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» siehe Seite 35

SICCO **Control Star, Control Big Star and Control Maxi 2 Desiccators**

OPERATING INSTRUCTIONS

1. SCOPE OF DELIVERY.....	14
2. SAFETY INSTRUCTIONS.....	14
3. INITIAL OPERATION.....	15
4. FUNCTIONAL NOTES.....	16
5. SERVICE AND MAINTENANCE.....	21
6. TROUBLESHOOTING / REPAIRS.....	21
7. DISPOSAL INSTRUCTIONS.....	21
8. SPECIFICATIONEN.....	22
9. SPARE PARTS / SPECIAL ACCESSORIES.....	23

SICCO **Control Star Desiccator**
Art. No.: V 1677-07

SICCO **Control Big Star Desiccator**
Art. No.: V 1679-07

SICCO **Control Maxi 2 Desiccator**
Art. No.: V 1682-07

Please read the following information carefully. Please strictly observe the safety, operating, and maintenance instructions.

SICCO Control-Star Desiccators are ideally suitable for long-term storage of reference samples, electronic components, photo equipment, papers, valuable cultural assets and much more. With the Control Star Desiccator, the humidity inside the desiccator can be kept permanently at an adjustable value (15 % to 95 % RH). The control unit uses the integrated silica gel to automatically maintain the target humidity within a tolerance of approx. $\pm 5\%$. The viewing panels are made of high-quality, transparent plastic; for the Vitrum model, they are made of borosilicate glass 3.3. The aluminium profiles are securely connected by polyamide connectors.

1. Scope of Delivery	Shelf	Guide Rails	Display	Operating Instructions
Control Star Desiccator	4 (made of PC/polycarbonate)	8 (made of PA-GF/polyamide with glass fibre)	1	1
Control Big-Star Desiccator	2 (made of stainless steel)	4 (made of PA/polyamide-black)	1	1
Control Maxi 2 Desiccator	4 (made of stainless steel)	8 (made of PA/polyamide-black)	2	1

2. Safety Instructions

- » Please do not operate a defective device.
- » SICCO Control Desiccators are only suitable for 100-230 V AC.
- » Please do not open the control box.
- » Please do not store solvents or acids in the desiccator.
- » Do not place explosive or flammable materials in the desiccator.
- » Do not generate overpressure or underpressure in the desiccator.
- » Ensure a safe and stable stand for the desiccator.
- » Only use the desiccator at a working height that is easily accessible to you.
- » Only use the desiccator at an ambient temperature of $-20\text{ }^{\circ}\text{C}$ to $+70\text{ }^{\circ}\text{C}$.
- » SICCO Control Desiccators are only suitable for operation in closed, dry rooms.
- » Protect the Control Desiccators from splashing water.
- » Protect the Control Desiccators from direct sunlight.
- » Pull out only one telescopic shelf at a time.
- » Please do not remove the rubber plug or the silica gel reservoir cover.

SICCO Control Star, Control Big Star and Control Maxi 2 Desiccators

3. Initial Operation

» Remove the desiccator from its protective cover.

» Remove the packaging from the shelves.

» To open the door, please hold the housing of the desiccator with one hand while opening the door by pulling the handle with the other hand.

» **Important information:** If a desiccator has been closed for a longer period of time, the force required to open the door may be greater than what you need in daily use. You can compensate for this by pulling the handle more strongly.


» Press the guide rails (1) horizontally and in pairs on both sides into the grid rails (2). The guide rail (1) audibly snaps into the grid rail (2). For better orientation, the grid rails (2) are numbered.


SICCO Control Star, Control Big Star and Control Maxi 2 Desiccators

- » Push a shelf (2) into the guide rails (1). Now mount another pair of guide rails (1) and push in another shelf (2).

Note: Assembling these in pairs and immediately inserting a shelf (2) prevents the guide rails (1) from accidentally coming loose.

- » Connect the small plug (1) to the socket (2) on the control unit. Plug the power supply (3) into a 100-230 V socket.


- » As soon as the control unit is supplied with power, the valve slides open and close once, the fan starts to work and the status LED lights up permanently green.


4. Functional Notes

A desiccator equipped with a control unit allows the targeted setting as well as the long-term maintenance of a defined humidity value in the interior.

- » **Dehumidification / Drying Operation:** In dehumidification mode, special desiccant pearls (e.g. silica gel) absorb the humidity of the interior (1). An integrated fan optimizes this process through a continuous air flow. For regeneration, a gate valve opens, the material is heated and the residual moisture is released to the outside ambient air.


- » **Humidification Operation:** In humidification mode (2), a fan guides ambient air past the desiccant pearls, which absorb its moisture. After a preset period of time, the gate valve closes. The desiccant is then heated and releases the humidity into the interior, whereby the relative humidity rises in a controlled manner.

The Display – Adjusting the Setpoint

The Control Star Desiccator display provides the following information:

- 1 Operating status (status)
- 2 Current temperature in °C or °F (temperature)
- 3 Actual value of relative humidity (actual)
- 4 Adjusted setpoint value of relative humidity (setpoint)
- 5 The setpoint value can be adjusted using the arrow keys
- 6 Select the temperature value to toggle between °C and °F.

1 Operating Status Indicators (status)

Display:	Function
Dehumidify:	The adjusted setpoint is lower than the current actual value. The control unit is in drying/dehumidification mode.
Humidify:	The adjusted setpoint is higher than the current actual value. The control unit operates in humidification mode.
OK:	The current humidity in the desiccator is within the setpoint tolerance range ($\pm 5\%$). The control unit is in monitoring mode; no humidification or dehumidification is occurring.

2 Temperature Display (temperature)

- » Shows the current temperature inside the desiccator in degrees Celsius or Fahrenheit.
- » By pressing on the temperature value, you can switch between °C and °F.

3 Actual Value („actual“)

- » Display of the currently measured humidity value in % relative humidity.

4 Setpoint Value („setpoint“)

- » Use the arrow keys to increase or decrease the humidity setpoint.
- » The values change in 1% steps.
- » Hold the key to adjust the value in 5% steps.

SICCO **Control Star, Control Big Star and Control Maxi 2 Desiccators**

The Control Unit

The coloured LED indicators signal the current operating status.


Left Unit (Viewed from the Front):

- » Responsible for reducing the humidity.
- » Permanently glowing green: Unit active, drying the interior.
- » Flashing green: Desiccant is regenerating.

Right Unit (Viewed from the Front):

- » Responsible for humidifying the interior.
- » Flashing blue: Humidification in progress.
- » Steady green: The unit is absorbing humidity from the environment.


LED off:

- » Setpoint reached – monitoring mode

Display


- » **Temperature Range:** -40°C to +85°C
- » **Humidity:** 15 to 95 %
- » **Temperature Accuracy:** +/- 1°C
- » **Humidity +/- 3 %**

Achievable Values: Min/Max

To determine the limits, setpoints of 0% and 100% relative humidity were used.
The humidity control trend was documented.


Minimum Values (Dehumidification)


Time in hours

Conclusion: After approx. 5 hours, the relative humidity already drops below 20 %.
After about 10 to 12 hours, a stable value of approx. 10 % is reached.

Maximum Values (Humidification)


Time in minutes

Time in minutes

Conclusion: After approx. 10 minutes, the relative humidity rises to over 90 %.
Subsequently, a humidity level exceeding 95% is reached.

Control Accuracy (Tolerance Range)

Time in hours

Conclusion: At a setpoint of 65 % relative humidity – outside the minimum and maximum range – a control accuracy of ± 5 % can be maintained over the entire test period.

5. Service and Maintenance

- » Please use only soft cloths for cleaning and do not use any solvent-based or abrasive cleaning agents.
- » Keep the magnetic strip of the door clean.

6. Troubleshooting / Repairs




The following troubleshooting instructions are intended to help you identify and resolve potential malfunctions. If you are unable to repair the damage yourself, please contact us by phone. We will be glad to help you resolve the issue. Should these steps fail to resolve the problem, please return the device cleaned and in its original packaging. We will then assess the unit to determine the most cost-effective repair options.

Contact our customer support for the following error messages!

Customer Support: +49 (0) 93 46-92 86-0

E-Mail-Support: info@sicco.de

Important: Please state model and serial number.

Complaint	Cause and Remedy
The humidity in the interior does not change.	The desiccator is not connected to the power supply. Plug the plug into a 100-230 V socket. The door is not closed tightly. Close the door. The built-in hygrometer is defective. Test the built-in hygrometer with another hygrometer.

7. Disposal Instructions

The desiccators and the packaging are made of recyclable materials and can therefore be disposed of at local waste disposal facilities. For more detailed information on the materials used, please refer to the following specifications.

8. Specificationen

SICCO Control Star Desiccators - Art. No.: V 1677-07

Drying performance:	approx. 70 g humidity in 24 hours
Humidification performance:	approx. 30 g humidity in 24 hours
Power consumption:	max. 120 W at 24 V
Mains voltage:	230 volts
Material of panels:	PMMA (polymethyl methacrylate), door and back panel made of PC (polycarbonate)
Temperature resistance:	20 °C to + 70 °C
Frame material:	aluminium with polyamide connectors
Locking mechanism:	Magnetic seal
Overall dimensions W x H x D:	310 x 525 x 375 mm
Drying unit extension:	60 mm at the back
Usable interior space W x H x D:	260 x 480 x 330 mm
Weight:	approx. 7,3 kg
Max. load per shelf/ total load:	10 kg / 30 kg

SICCO Control Big Star Desiccators - Art. No.: V 1679-07

Drying performance:	approx. 115 g humidity in 24 hours
Humidification performance:	approx. 80 g humidity in 24 hours
Power consumption:	max. 240 W at 24 V
Mains voltage:	230 volts
Panel material:	PMMA (polymethyl methacrylate), door and back panel made of PC (polycarbonate)
Temperature resistance:	20 °C to + 70 °C
Frame material:	aluminium with polyamide connectors
Locking mechanism:	Magnetic seal
Overall dimensions W x H x D:	560 x 560 x 580 mm
Drying unit extension:	60 mm at the back
Usable interior space W x H x D:	495 x 500 x 540 mm
Weight:	approx. 23 kg
Max. load per shelf/ total load:	10 kg / 80 kg

SICCO Control Maxi 2 Desiccators - Art. No.: V 1682-07

Drying performance:	approx. 115 g humidity in 24 hours per chamber
Humidification performance:	approx. 80 g humidity in 24 hours per chamber
Power consumption:	max 2x 240 W at 24 V
Mains voltage:	230 volts
Material of panels:	PMMA (polymethyl methacrylate), door and back panel made of PC (polycarbonate)
Temperature resistance:	20 °C to + 70 °C
Frame material:	aluminium with polyamide connectors
Locking mechanism:	Magnetic seal
Overall dimensions W x H x D:	560 x 1150 x 580 mm
Drying unit extension:	60 mm at the back
Usable interior space W x H x D:	495 x 500 x 540 mm per chamber
Weight:	approx. 36 kg
Max. load per shelf/ total load:	10 kg / 60 kg

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For the following products:

Device type: Drying Cabinets

Type designation: Control Star Desiccator	Cat. No. V 1677-07
Control Big Star Desiccator	Cat. No. V 1679-07
Control Maxi 2 Desiccator	Cat. No. V 1682-07

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Hereby declares that the aforementioned products comply with the requirements of the following directives:

2014/30/EU EMC Directive (Electromagnetic Compatibility)

2014/35/EU Low Voltage Directive

The products are **RoHs** compliant. They comply with Directive

2011/65/EU Directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment, including **amendment 2015/863**.

To ensure compliance, the following standards/specifications were referenced:

EN 61000-6-2 2019

EN 61000-6-4+A1 2011

EN 62368-1 2014

EN 61326-1 2012

EN 61010-1 2019

EN IEC 63000 2019-05

The product fulfills the requirements provided that its intended use is ensured during installation and operation.

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
BOLA

Special Requirements? Custom-Made

You are looking for something very particular? Somthing that even our giant portfolio of thought-through solutions does not cover?

No problem. As developer and producer we offer you the pdsibility of an individual production, as desired. This can be made faster and more economical than expected.

Just talk to our specialists about your ideas - we advise and support you already beginning with the construction and produce than material-appropriate according to your specifications. And this already from numer of pieces 1 at an inexpected price than you imagine.

For this, we only need a drawing (a rough sketch is sufficient) and some information.

The check list for your tailor-made dessiccator/ glove box:


- ✓ What is the name of your dessiccator/ glove box?
- ✓ For which application do you need the desiccator/ glove box?
- ✓ Which outer dimensions do you need?
- ✓ What do you want to store?
- ✓ Which panel material do you require?
- ✓ What kind of shelves do you require?
- ✓ What will be the maximum weight load per shelf?
- ✓ What kind of accessories do you require (cylinder look, cable lead-ins, etc.)?
- ✓ Which number of pieces do you require?
- ✓ What is the budgetary price range per desiccator that should not be exceeded?

Do you have a special request?
Just call us: **+49 (0) 93 46 92 86-0**

Or just send us a drawing (rough scetch is sufficient) and some information to **info@sicco.de**. We will come back to you to discuss the details and establish a free offer.

Desiccators and Drying Cabinets

Drying Cabinets, Desiccators and Glove Boxes

Stable | Functional | Custom-Made

For work in an isolated atmosphere and for dust- and moisture-free storage or storage with inert gas.

Ask for your catalogue or get more information:

www.sicco.de/en or +49 (0) 93 46-92 86-0

A PRODUCT BRAND OF BOHLENDER

sicco.de/en

You can find more of our
products on sicco.de/en

Desiccators

Gloves

Desiccant

Drying Cabinets

Hygrometer

Shelves

SICCO

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BOHLENDER GmbH

Waltersberg 8

D 97947 Grünsfeld

+49 (0) 93 46-92 86-0

info@sicco.de

www.sicco.de