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Bedienungsanleitung Operating Instructions Mode d'emploi

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
BOLA **GL-Kugelhähne**

OPERATING INSTRUCTIONS

1. SCOPE OF DELIVERY	8
2. SECURITY INSTRUCTIONS	8
3. MODE OF FUNCTIONS	8
4. TAKING INTO OPERATION	8
5. SERVICE / MAINTENANCE	9
6. DISTURBANCES / REPAIRS	12
7. WASTE DISPOSAL	13
8. SPECIFICATIONS	13
9. SPARE PARTS / ACCESSORIES	13

BOLA **GL Ball Valve,** **Type 2-Way**

Cat.No.:

- E 664-10
- E 664-20
- E 664-30
- E 664-40

BOLA **GL Ball Valve,** **Type 3-Way**

Cat.No.:

- E 667-10
- E 667-20
- E 667-30
- E 667-40

BOLA GL Ball Valves

Please read the information listed below thoroughly. The given information about security, use and maintenance has to be observed. BOLA Ball Valves can be used with virtually all kinds of media that do not contain any solids. They are suitable for vacuum and can be used under pressure up to 12 bar .

1. Scope of Delivery

- 1 Ball Valve
- 1 Operating Instruction

2. Security Instructions

- » Do not take a damaged valve into operation.
- » Use the ball valves only for gas or watery mediums that do not contain any solids.
- » The working pressure at room temperature has to be limited to max. 12 bar.
- » Do only apply working temperatures from -20 to +220°C.
- » Before connecting or disconnecting the valve, please make sure that there is no pressure/vacuum on the supply pipes.

3. Mode of Function

By turning the handle, the ball stopcock plug can be moved to the requested position and the valve be opened or closed. The flow direction of the 3-way ball valve is defined by an arrow on the handle. Connect the inlet to the lower socket. The 2-way ball valve is open when the handle is in parallel position to the thread sockets.

4. Taking into operation

- » Integrate the ball valve with a tubing to your line system
- » Use an appropriate laboratory screw joint (20), for example from BOLA (not included in delivery)
- » 2-Way Ball Valve: you don't have to pay attention on the flow direction

- » Integrate the ball valve with a tubing to your line system
- » Use an appropriate laboratory screw joint (20), for example from BOLA (not included in delivery).
- » 3-Way Ball Valve: connect your supply pipe to the lower socket (30).

BOLA GL Ball Valves

5. Service / Maintenance


Do not use any abrasive cleaning agents or sharp tools for cleaning. Rinse the valve after use with clear water. Open and close the valve several times during rinsing. The valve can be disassembled for cleaning as follows:

- » Screw off the sealing bolts in the threaded sockets with a screwdriver by turning them counter-clockwise.
- » If a 3-way ball valve (E 667-..) is disassembled, only the lateral sealing bolts have to be removed.


Attention: Please take care not to glide off – risk of injury!

- » Screw off the black screw cap (2) by turning it counter-clockwise and remove it completely from the thread. Pull out carefully the ball stopcock plug assembly (4) of the base body (3).


- » **Attention! Do not damage any sealing surfaces (marked with arrows) . If sealing surfaces are damaged, the functionality of the ball valve cannot be guaranteed anymore!**

All single components can be easily cleaned with usual, non-abrasive cleansing agents. The assembly after cleaning is carried out vice versa. **Please pay special attention to the following instructions:**

BOLA GL Ball Valves

- » Put the ball stopcock plug assembly (4) with slight pressure into the base body (3).


- » Screw the black screw cap (2) one time counter-clockwise on the thread of the base body.

- » Tighten the screw cap (2) by turning it clockwise until the passage of the stopcock plug is parallel to the passage of the base body.

- » Put the sealing bolts (1) into the sockets (3).

- » Tighten the sealing bolts with a screw driver by turning them clockwise until you notice a light drag.

Attention: Please take care not to glide off – risk of injury!!!


BOLA GL Ball Valves


All single components can be easily cleaned with usual, non-abrasive cleansing agents. The assembly after cleaning is carried out vice versa.

Please pay special attention to the following instructions:

» Put the ball stopcock plug assembly (4) with slight pressure into the base body (3).


» Screw the black screw cap (2) one time counter-clockwise on the thread of the base body.

» Tighten the screw cap (2) by turning it clockwise until the passage of the stopcock plug is parallel to the passage of the base body.


» Put the sealing bolts (1) into the sockets (3).

» Tighten the sealing bolts with a screw driver by turning them clockwise until you notice a light drag.

Attention:

Please take care not to glide off – risk of injury!!!

BOLA GL Ball Valves

» After assembly, the sealing bolts have to be adjusted.

All sealing bolts have to be screwed in at equable depth. Verify the depths by measuring with a calliper the depth from the sealing surface of the socket to the sealing bolt. The screw-in depths shall not deviate more than 0,2 mm.

If the ball stopcock plug is stuck or can hardly be moved, the sealing bolts (1) have to be equably unscrewed. Check the screw-in depths with a calliper after every new adjustment.

» Tighten the screw cap (2) manually without any tools by turning it clockwise.

6. Disturbances / Repairs

The following information shall help you to recognise and, if possible, to eliminate possible reasons for disturbances by yourself. If you are not able to handle the disturbance, please give us a call. We will then try to help you this way. If then, a solution to your problem cannot be found, we will ask you to return the cleaned BOLA Ball Valve to us. Then we will check whether a cost-effective repair is possible.

Problem	Reason and Help
No flow:	<ul style="list-style-type: none"> » Seed-line is blocked. Remove the blockage or use a new feed-line. Passage of stopcock plug is blocked. Disassemble the ball valve and remove the ball stopcock plug assembly. You can now clean thoroughly the stopcock plug. » The valve is closed. Please open it.
Hahn ist undicht:	<ul style="list-style-type: none"> » Laboratory screw joint is not mounted correctly. Disassemble the laboratory screw joint and check whether the inner parts are mounted in the correct order (see below). <div style="text-align: center; border: 1px solid red; padding: 5px; margin: 10px 0;"> </div> <ul style="list-style-type: none"> » The sealing bolts are worn. Please order new sealing bolts (see 5). » The sealing surface on the GL-socket has been damaged. Please return the cleaned valve to the manufacturer for inspection. » The screw cap (2) of the ball valve is not tightened. Tight the screw cap by turning it clockwise. Please take care not to damage the thread by tighten it too strong.
The valve can neither be opened nor closed:	<ul style="list-style-type: none"> » Tubing might be pulled too deep into the valve and blocks the stopcock plug. Disconnect the laboratory screw joint and pull the tubing out of the stopcock plug to remove the blockage.


BOLA GL Ball Valves

7. Waste Disposal


The valve as well as the packing consist of recyclable materials.

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BOHLENDER GmbH