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
Bedienungsanleitung Operating Instructions Mode d'emploi

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BOLA **Non-return valves**

OPERATING INSTRUCTIONS

1. SUPPLIED WITH	10
2. SECURITY INSTRUCTIONS	10
3. MODE OF FUNCTION	10
4. TAKING INTO OPERATION	11
5. SERVICE / MAINTENANCE	13
6. DISTURBANCES/ REPAIRS	14
7. WASTE DISPOSAL	15
8. SPECIFICATIONS	15
9. SPARE PARTS / ACCESSORIES	15

BOLA Non-return valves
Cat-No.: E 680-21, -23

BOLA Non-return valves
Cat-No.: E 680-27, -31

Please read the information listed below thoroughly. We kindly ask you to observe the information about security, use and maintenance.

BOLA Non-Return Valves can be connected to tubing or tubes. The integrated lock function prevents a backflow of the liquid. The opening pressure in flow direction can be adjusted between 0.1 and 2.0 bar. The flowing liquid or gas is only exposed to PTFE and PFA.

BOLA **Non-return valves**

1. Supplied with

- » 1 BOLA Non-return valve
- » 2 Screw nuts
- » 2 Pairs of V-rings and sealing rings
- » 1 Operating Instructions

2. Security instructions

- » Please do not take a damaged valve into operation.
- » Only use the valve for substances in a temperature range between -200°C and +250°C.
- » Only use the valve within a pressure range of 0.1 and 2.0 bar.
- » Only use the valve for gases or watery substances without solids.
- » Please consider the flow direction. A corresponding arrow is engraved on the housing.
- » Make sure that all components of your system resist the pressure that you have adjusted in the valve.
- » Check the requested opening pressure before integrating the valve in your system.
- » Inlet and outlet tubing must be pressureless before integrating the valve in or removing it from your system.

The operating instruction cannot show all security problems that can occur. It is the user's responsibility to observe security and health instructions and to fix corresponding restrictions before use.

3. Mode of function


A free flow is only possible in flow direction (engraved arrow on housing). The opening pressure is adjusted by a setscrew inside the valve. In flow direction, the piston is lifted by the liquid/gas and an annular passage allows the liquid/gas to flow. This passage is closed as soon as the pressure falls below the adjusted value. The closing pressure is caused by a spring so that the fitting position is arbitrary.

BOLA **Non-return valves**


4. Taking into operation

The valve is supplied with an opening pressure of 0.1 bar. If you want to adjust another pressure value, please proceed as follows:


- » Remove the nut (1) on the outlet side by turning it in counterclockwise direction and remove the V-ring and sealing ring. The engraved arrow points to the outlet side.
- » Open the nut (2) on the inlet side by making approx. two turns in counterclockwise direction.

- » Connect the valve with a tubing (20) on the inlet side to an air service unit. Push the tubing to the bottom of the screw joint.
- » Tighten the nut of the screw joint (2) by turning in clockwise direction with two fingers.

- » Please consider: The O.D. of a tubing can vary due to production tolerance. If the tubing cannot be pushed into the screw joint, we recommend to assemble the components separately. Please observe the arrangement of sealing ring (3) and tapered ring (4).


BOLA Non-return valves are available for different tubing diameters. For the following adjustments, different sizes of screwdrivers are needed:

- » For non-return valves for tubing O.D. 4 mm and 6 mm use a screwdriver with a 3.5 mm tip.
- » For non-return valves for tubing O.D. 8 mm and 10 mm use a screwdriver with a 5.5 mm tip.
- » For non-return valves for tubing O.D. 12 mm use a screwdriver with a 10 mm tip.

BOLA **Non-return valves**


- » Position a suitable screwdriver in the groove of the setscrew.
- » Turn the setscrew in clockwise direction as far as it will go. The non-return valve is now adjusted to the highest possible opening pressure.

- » Charge the non-return valve with the requested pressure (must be between 0.1 and 2.0 bar).

- » Position the suitable screwdriver into the groove of the setscrew and immerse it in a vessel filled with water.
- » Slowly turn the setscrew in counterclockwise direction until air bubbles come out of the valve. The opening pressure of the valve is now adjusted to the value of your air service unit.


- » The engraved arrow shows the flow direction of the non-return valve. Please consider this when integrating the valve into your system.


5. Service / Maintenance

Rinse the valve after use with a suitable cleaning solution (for example warm water). Connect the valve in flow direction for cleaning. The adjusted opening pressure possibly has to be reduced since it must be lower than the pressure of the inlet pipe. If you want to disassemble the valve for cleaning, please proceed as follows:


- » Remove the connecting nut (3) from the housing (4) by turning in counterclockwise direction.

- » Remove the connecting piece (5), the piston (6) and the spring (7).

- » Unscrew the setscrew with a suitable screwdriver completely from the housing (4) by turning in counterclockwise direction.


- » We recommend not to disassemble the o-ring inside the piston (6).


BOLA **Non-return valves**

All components can be cleaned with common, non-abrasive detergents. Assembly must be made in reverse order. Please observe following remarks:


- » Screw the setscrew (8) with a suitable screwdriver in the housing (4) by turning in clockwise direction. The counter bore of the setscrew (8) points outwards.


- » Push the connecting piece (5) into the connecting nut (3).

- » Position the piston (6) into the spring (7). Push the piston (6) and the spring (7) into the housing (4). Centre the spring (7) in the counterbore of the setscrew (8).

- » Screw the connecting nut (3) on the housing (4) by turning in clockwise direction. Tighten the connecting nut by hand. The requested opening pressure can be adjusted as mentioned above.

6. Disturbances/ Repairs

The following information shall help you to find by yourself the reason and a possible solution for disturbances. In case of disturbances you cannot handle, please do not hesitate to contact us by phone. Of course we will try to help you and give advice to find a solution. If the problem might not be solved this way, the cleaned valve can be shipped to us. We will check then if a cost-effective repair is possible.

BOLA **Non-return valves**

Problem	Reason and help
Liquid/gas does not flow.	» The flow direction was not observed. Check connections and change if necessary. The flow direction is marked with an arrow on the housing. » The adjusted opening pressure of the non-return valve is too high. It must be lower than the pressure of your system. Reduce the opening pressure of your non-return valve.
Liquid/gas flows contrary to flow direction.	» The valve is blocked by solids. Prevent this by assembling a filter. » The liquid is too viscous. The non-return valve is not suitable for viscous liquids.

7. Waste disposal

The BOLA Non-return valve as well as the packing consist of recyclable materials. Please take further information about the materials used from the following specifications.

8. Specifications

BOLA Non-return valve

Cat. No.	Thread of fitting	For tubing O.D.	Materials	Temperature resistance	Opening pressure	Adjusted pressure
E 680-21	M 14 x 2	4 mm	PTFE Glass-filled PTFE PFA Perfluor-elastomer Spring steel 1.4310	- 200 °C to + 250 °C	0,1 to 2,0 bar	0,1 bar
E 680-23	M 14 x 2	6 mm				
E 680-27	M 18 x 2	8 mm				
E 680-31	M 18 x 2	10 mm				

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BOHLENDER GmbH