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» siehe Seite 30-31

OPERATING INSTRUCTIONS

1. SUPPLIED WITH.....	12
2. SECURITY ADVICE.....	12
3. TAKING INTO OPERATION.....	12
4. FUNCTIONALITY.....	14
5. SERVICE / MAINTENANCE.....	15
6. DISTURBANCES / REPAIRS.....	15
7. WASTE DISPOSAL.....	16
8. SPECIFICATIONS.....	16
9. ACCESSORIES.....	30

SICCO Star-Desiccator
Cat.-No.: V 1871-07

SICCO Tower-Star-Desiccator
Cat.-No.: V 1938-07

SICCO Star-Desiccator Horizontal
Cat.-No.: V 1899-07

SICCO Star-Desiccator Antistatik
Cat.-No.: V 1910-07

SICCO Star-Desiccator Protect
Cat.-No.: V 1879-07

SICCO Star-Desiccator Black
Cat.-No.: V 1891-07

SICCO Star-Vitrum-Desiccator
Cat.-No.: V 1841-07

SICCO Star-Desiccator-Double
Cat.-No.: V 1930-07

SICCO Star-Vitrum-Desiccator Horizontal
Cat.-No.: V 1821-07

Please read the information listed below thoroughly. We urgently ask you to observe the instructions about security, use and maintenance.

SICCO Desiccators are made of high-quality and transparent plastic material, except the models VITRUM with panels made of borosilicate glass 3.3. They are made for safe and secure storage of humidity-sensitive products as well as of optical and microelectronic parts in a dust-free atmosphere. This happens by adding Silicagel which considerably reduces the air humidity by its absorbing effect.

1. Supplied with

- 1 SICCO Star-Desiccator
- 4 Shelves (Versions Horizontal: 2 shelves)
- 8 Rails for shelves (Versions Horizontal: 4 rails)
- 1 Hygrometer (fitted into the door)
- 1 Desiccant tray
- 1 Bag with 500 ml Silicagel (orange)
- 1 Operating Instructions

2. Security Advice

- » Do not take a damaged desiccator into operation.
- » Do not store any solvents or acids in the desiccator except for Star Vitrum Desiccators.
- » Do not put any explosive or inflammable products into the desiccator.
- » SICCO Desiccators are not suitable for overpressure or vacuum.
- » Provide a safe stand for the desiccator.
- » Place the desiccator in a reachable working height.
- » The ambient temperature range of $-20\text{ }^{\circ}\text{C}$ to $+70\text{ }^{\circ}\text{C}$ respectively $-70\text{ }^{\circ}\text{C}$ to $+150\text{ }^{\circ}\text{C}$ for models Vitrum shall not be exceeded.
- » Do not expose the desiccator to direct sunlight.
- » Remove only one shelf at once.

3. Taking into operation

- » Remove the protection cover.

- » For models with aluminium rails: Remove the protection film from the rails.

SICCO Star-Desiccators

- » Remove the protection film from the shelves.

- » Opening the door: Hold the desiccator housing with one hand while you open the door by pulling the door grip with the other hand. Please note: Due to the strong magnetic lock, it is harder to open the door the longer it was closed before. This is normal and can be compensated by pulling stronger than in the everyday use.

- » For models with aluminium rails: Hook the rails (2) into the lateral mounted suspension arrangements (snap-in pin (1) upwards).

- » For models with aluminium rails: Mount the rails (2) horizontally in pairs. Pay attention to mount them in the same height. A slight rotation of the rail helps to fix the snap-in pin safely in the suspension arrangement (3).


- » For models with rail (2) and suspension arrangement (3) made of plastic: Mount the rails (2) horizontally in pairs on both sides of the suspension arrangement (3). The engagement of the rail (2) is clearly audible. For a better orientation, the levels are marked by numbers on the suspension arrangement (3).

- » The shelves (4) can then be slid on the rails (2).
- » Mount another pair of rails (2) and slide in another shelf (4).

Please note: Mount the rails in pairs and slide in immediately a shelf to avoid that other rails can fall down accidentally.

SICCO Star-Desiccators

- » Open the bag with the Silicagel.

- » Fill the tray (included in delivery) with Silicagel. The desiccant can be used immediately.

- » Open the battery case on the back of the hygrometer by pushing the cover downwards.
- » **Attention (only for Star-Desiccator Black):** The Hygrometer cannot be exchanged.
- » Remove the battery safety strip and reinsert the battery. Close the battery case.

The following functions are only available with SICCO Star-Desiccator Antistatic Cat.-No. V 1910-07

- » Connect the ground cable on the back of the desiccator.

4. Functionality

The included Silicagel is very hygroscopic and absorbs all air humidity inside the desiccator. The Silicagel contains a colour indicator. In unsaturated state, the Silicagel is orange. With increasing saturation, the colour is changing from dark brown to black.

The Silicagel can be regenerated by putting it into an oven for 90 minutes at about 90 to 110 °C. The dumping height influences the time for regeneration: each centimetre in height will take about 90 minutes for regeneration. A microwave oven is not suitable!

Functions of the digital Thermometer/Hygrometer

- Measuring range:
- » Temperature – 10 °C to + 60 °C.
 - » Air humidity 10 % to 99 % relative humidity.
- Max. error in measurement:
- » +/- 1°C at –10 °C to + 50 °C.
 - » +/- 5 % relative humidity at 25 % to 75 % relative humidity and 0 °C to + 50 °C.
 - » **Attention:** The mounted thermometer/hygrometer cannot display any values or data if the relative inside humidity of the desiccator is less than 25 % or higher than 75 %. To display exact values, we recommend using a thermometer / hygrometer of higher quality.
- Battery:
- » Please use a button cell LR44 or AG13.
- Functionality of the Keys:
Indication of maximum / minimum values:
- » By pressing the MAX/MIN-key, the maximum value since the last reset will be displayed.
 - » By pressing the MAX/MIN-key again, the minimum value since the last reset will be displayed.
 - » The current values will be displayed after pressing MAX/MIN a third time.
 - » Reset of the MAX/MIN-memory:
Press CLEAR and the MAX/MIN-memory will be deleted. The current values will be stored instead.
 - » Battery sign in the display
If the battery is low, a battery flag is shown on the left bottom of the display. Exchange the battery to ensure best possible operation of the device.

5. Service / Maintenance

- » Use only soft clothes for cleaning. Do not use any detergents which may contain solvents or which scrub.
- » Keep the magnetic door strip clean.

6. Disturbances / Repairs

The following information shall help you to recognise and, if possible, to eliminate possible reasons for disturbances by yourself. If you are not able to handle the disturbance, please give us a call. We will then try to help you this way. If then, a solution to your problem cannot be found, we will ask you to return the cleaned desiccator in the original packing to us. Then we will check whether a cost-effective repair is possible.

Objection

Humidity inside the desiccator does not fall.

Reason and help

- » The door is not closed tightly. Close the door.
- » The battery in the hygrometer is empty. Exchange the battery.
- » The Silicagel is saturated. Regenerate the Silicagel as described under point 4.
- » Perhaps you have not added enough Silicagel. Add further Silicagel.
- » The built-in hygrometer might be faulty. Double-check the displayed values with an external hygrometer. Exchange the hygrometer.

7. Waste Disposal

The desiccators as well as their packing are made of recyclable materials. Please take further information on the materials used from the following specifications.

8. Specifications

SICCO Star-Desiccator

Cat.-No:	V 1871-07
Material panels:	PMMA (Polymethylmethacrylate)
Temperature resistance:	- 20° to + 70°C
Material frame:	Aluminium with connectors made of polyamide
Material of rails:	Polyamide
Lock:	Magnetic tape
Overall dimensions (W x H x D):	310 x 525 x 375 mm
Usable interior space (W x H x D):	260 x 480 x 330 mm
Weight:	approx. 7 kg
Max. all-over load per shelf:	10 kg
Total all-over load:	30 kg

SICCO Tower-Star-Desiccator

Cat.-No.:	V 1938-07
Material panels:	PMMA (Polymethylmethacrylate)
Temperature resistance:	- 20° to + 70°C
Material frame:	Aluminium with connectors made of polyamide
Material of rails:	Aluminium
Lock:	Magnetic tape
Overall dimensions (W x H x D):	310 x 735 x 375 mm
Usable interior space (W x H x D):	260 x 690 x 330 mm
Weight:	approx. 8 kg
Max. number of shelves:	38 pieces
Max. all-over load per shelf:	10 kg
Total all-over load:	40 kg

SICCO Star-Desiccator Horizontal

Cat.-No:	V 1899-07
Material panels:	PMMA (Polymethylmethacrylate)
Temperature resistance:	- 20° to + 70°C
Material frame:	Aluminium with connectors made of polyamide
Material of rails:	Polyamide
Lock:	Magnetic tape
Overall dimensions (W x H x D):	525 x 340 x 375 mm
Usable interior space (W x H x D):	480 x 260 x 330 mm
Weight:	approx. 7 kg
Max. number of shelves:	13 pieces
Max. all-over load per shelf:	1.5 kg
Total all-over load:	30 kg

SICCO Star-Desiccator Antistatic

Cat.-No:	V 1910-07
Material panels:	PC-ESD
Temperature resistance:	- 20° to + 70°C
Material frame:	Aluminium with connectors made of polyamide
Material of rails:	Aluminium
Lock:	Magnetic tape
Overall dimensions (W x H x D):	310 x 525 x 375 mm
Usable interior space (W x H x D):	260 x 480 x 330 mm
Weight:	approx. 7 kg
Max. number of shelves:	26 pieces
Max. all-over load per shelf:	10 kg
Total all-over load:	30 kg

SICCO Star-Desiccator Protect

Cat.-No.:	V 1879-07
Material panels:	PMMA orange
Temperature resistance:	- 20° to + 70°C
Material frame:	Aluminium with connectors made of polyamide
Material of rails:	Polyamide
Lock:	Magnetic tape
Overall dimensions (W x H x D):	310 x 525 x 375 mm
Usable interior space (W x H x D):	260 x 480 x 330 mm
Weight:	approx. 7 kg
Max. number of shelves:	26 pieces
Max. all-over load per shelf:	10 kg
Total all-over load:	30 kg

SICCO Star-Desiccator Black

Cat.-No:	V 1891-07
Material panels:	PMMA black
Temperature resistance:	- 20° to + 70°C
Material frame:	Aluminium with connectors made of polyamide
Material of rails:	Polyamide
Lock:	Magnetic tape
Overall dimensions (W x H x D):	310 x 525 x 375 mm
Usable interior space (W x H x D):	260 x 480 x 330 mm
Weight:	approx. 7 kg
Max. all-over load per shelf:	10 kg
Total all-over load:	30 kg

SICCO Star-Vitrum-Desiccator

Cat.-No.:	V 1841-07
Material panels:	Borosilicate glass 3.3
Temperature resistance:	- 70° to + 150°C
Material frame:	Aluminium with connectors made of polyamide
Material of rails:	Aluminium
Lock:	Magnetic tape
Overall dimensions (W x H x D):	310 x 525 x 375 mm
Usable interior space (W x H x D):	260 x 480 x 330 mm
Weight:	approx. 14 kg
Max. number of shelves:	26 pieces
Max. all-over load per shelf:	10 kg
Total all-over load:	30 kg

SICCO Star-Vitrum-Desiccators Horizontal

Cat.-No:	V 1821-07
Material panels:	Borosilicate glass 3.3
Temperature resistance:	- 70° to + 150°C
Material frame:	Aluminium with connectors made of polyamide
Material of rails:	Aluminium
Lock:	Magnetic tape
Overall dimensions (W x H x D):	525 x 340 x 375 mm
Usable interior space (W x H x D):	480 x 260 x 330 mm
Weight:	approx. 12 kg
Max. number of shelves:	13 pieces
Max. all-over load per shelf:	10 kg
Total all-over load:	30 kg

SICCO Star-Desiccator Double

Cat. No.:	V 1930-07
Material panels:	PMMA (Polymethylmethacrylate)
Temperature resistance:	- 20° to + 70°C
Material frame:	Aluminium with connectors made of polyamide
Material of rails:	Polyamide
Lock:	Magnetic tape
Overall dimensions (W x H x D):	630 x 520 x 375 mm
Usable interior space (W x H x D):	570 x 480 x 330 mm
Weight :	approx. 13 kg
Max. number of shelves:	26 pieces
Max. all-over load per shelf:	10 kg
Total all-over load:	30 kg

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BOHLENDER GmbH