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Bedienungsanleitung Operating Instructions Mode d'emploi

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BOLA Temperature Probes PT 100

OPERATING INSTRUCTIONS

1. SUPPLIED WITH	11
2. SECURITY INSTRUCTIONS	11
3. TAKING INTO OPERATION	12
4. MODE OF FUNCTION	13
5. SERVICE / MAINTENANCE	14
6. DISTURBANCES / REPAIRS	14
7. WASTE DISPOSAL	14
8. SPECIFICATIONS	15
9. SPARE PARTS / ACCESSORIES	26
10. CUSTOMIZED TEMPERATURE PROBES	17
CE DECLARATION OF CONFORMITY	29

BOLA Temperature Probes
LEMO® SOCKET
PT 100 Compact
Cat.-No.: P 1730-...

BOLA Double-Temperature Probes
LEMO® SOCKET
PT 100 Compact
Cat.-No.: P 1740-...

BOLA Temperature Probes
KABEL STRANDS
PT 100
Cat.-No.:
P 1750-...

BOLA Temperature Probes
LEMO® SOCKET
PT 100
Cat.-No.: P 1760-...

BOLA Temperature Probes
LEMO® PLUG
PT 100
Cat.-No.: P 1762-...

BOLA Immersion Probes
KABEL STRANDS
PT 100
Cat.-No.: P 1780-20

BOLA Immersion Probes
LEMO® SOCKET
PT 100
Cat.-No.: P 1790-20

BOLA Immersion Probes
LEMO® PLUG
PT 100
Cat.-No.: P 1792-20

Please read the information listed below thoroughly.

We urgently ask you to observe the information about security and use.

BOLA Temperature Probes are PT 100 probes. They are ideal for measuring temperatures in aggressive liquids. They are also available with connectors type LEMO® either plug or socket.

The "LEMO® compact" probe as well as the "Double probe LEMO® compact" have one respective two LEMO® sockets, which are mounted directly on the sensor. The LEMO® connectors are already mounted and can be connected to every suitable measuring instrument.

1. Supplied with

1 Temperature Probe PT 100 with cable length 1.500 mm or Immersion Probe with cable length 4.000mm (both alternatively with or without LEMO® socket size 1 / LEMO® plug size 1) respectively Temperature Probe LEMO® Compact (single or double) with LEMO® socket size 1

2. Security instructions

The temperature probes are designed to make measurements in aggressive liquids. The operating instructions cannot show all security problems which can arise eventually. It is up to the user to make sure that all safety and health regulations are kept and that the corresponding restrictions are stipulated before use.

- » Please do not take a damaged temperature probe into operation.
- » Please take the response times of the temperature probe into consideration during use. Typical response times T 50 and T 90 (50% respectively 90% of the final measuring value) are to be found under point 8 Specifications on pages 15-17.
- » Please do not exceed the operating temperature range of maximum +250 °C and minimum -50 °C.
- » Please verify that the LEMO® connectors are linked correctly to the measuring instrument.

3. Taking into operation


- » The temperature / immersion probes only have to be connected to a suitable measuring instrument. After that, you can start measuring.


- » For the LEMO® compact probes (single and double) please use an adequate cable to connect the probe to the measuring instrument.

- » Temperature probes and immersion probes with cable strands allow to connect an alternative connector.

- » For temperature and immersion probes without LEMO® connectors, the following connecting plan has to be considered. If you want to mount your own connector, please consider the correct position of the pins.

4. Mode of function

BOLA temperature probes are sensors based on platinum resistance changes under temperature influence. They are manufactured in conformance to class A. Thus, depending on the real temperature T, the allowed temperature tolerance dT is as follows: the deviation of the measured temperature is equal or less than $\pm 0,15 \text{ }^\circ\text{C} + (0,002 \times T)$.

The sensor itself is located at the end of a PTFE-encapsulated stainless steel tube (material code 1.4571). This tube provides certain rigidity, but can manually be bent to the requested form. The seamless PTFE-encapsulation provides an almost universal chemical resistance.

» BOLA Temperature Probe PT 100

The probes have tapered tips which reduce the response time considerably. The cable is encapsulated with PFA and connected tightly to the temperature probe. The 4-wire cable has a length of 1.500 mm and is equipped with a LEMO® connector, optionally available as socket size 1 or plug size 1.

» BOLA Immersion Probe PT 100


This probe can be immersed completely in aggressive liquids and gases. The 4-wire cable is coated with white PFA. The cable and the probe are resistant to virtually all chemicals. In addition, the total immersion probe can be used for measurement of temperature in autoclaves. They are equipped with a LEMO® connector, optionally available as socket size 1 or plug size 1 and 4-wire cable with a length of 4.000 mm.

» BOLA Temperature Probe LEMO® Compact PT 100

The probes have tapered tips which reduce the response time considerably. To connect a cable in any requested length, the LEMO®-socket size 1 has been fixed directly at the end of the probe.

» BOLA Double Temperature Probe LEMO® Compact PT 100

This probe is equipped with two independent thermo elements (1) and (2). To connect cables in any requested lengths, the LEMO®-sockets size 1 have been fixed directly at the end of the probe.

5. Service / Maintenance

In case of a strong contamination, the temperature probe should be cleaned after use. Generally, all fluoro-plastic materials, PTFE, PFA and FEP, have smooth, non-wetting surfaces and can usually be cleaned without any problems. Please do not use any abrasive scouring agents as they might damage the surfaces. You may use all neutral detergents (pH 7). For a stronger contamination, we recommend the use of an alkaline detergent up to pH 12. Apart from that, the probes are maintenance-free.

6. Disturbances / Repairs

If the probe should not work properly, please check first if the LEMO® connector is linked correctly to the measuring instrument. Also a damage of the connecting cable or of the sensor can lead to disturbances. Should there be a disturbance you cannot handle, please do not hesitate to call us. Of course we will try to help you this way. If we do not find a solution to your problem then, we will ask you to return the cleaned probe in the original packing to us. We will check if a cost-effective repair is possible.

7. Waste disposal

The probes as well as their packing are made of recyclable materials. Please take further information about the materials used from the following specifications.

8. Specifications » BOLA Temperature Probes PT 100

BOLA Cat. No.	P 1750-...
Materials	Stainless steel tube 1.4571 encapsulated with PTFE, PFA-coated cable
Temperature range	max. +250°C, min. -50°C
Resistance	100 Ohm at 0 °C
Probe-Ø	8 mm
Useful length	100 / 200 / 300 / 500 / 600 mm
Total length, approx	160 / 260 / 360 / 560 / 660 mm
Tip-Ø	5,5 mm
Collar ring	approx. 12 mm
Cable length	1.500 mm
Response time	T 50: 7-12 s; T 90: 14-16 s

» BOLA Temperature Probes PT 100 LEMO® either with socket size 1 or plug size 1

BOLA Cat. No.	P 1760-... / P 1762-...
Materials	Stainless steel tube 1.4571 encapsulated with PTFE, PFA-coated cable, Chromed brass (LEMO® coupling)
Temperature range	max. +250°C, min. -50°C
Resistance	100 Ohm at 0 °C
Probe-Ø	8 mm
Useful length	100 / 200 / 300 / 500 / 600 mm
Total length, approx.	160 / 260 / 360 / 560 / 660 mm
Tip-Ø	5,5 mm
Collar ring-Ø	approx. 12 mm
Cable length	1.500 mm
Response time	T 50: 7-12 s; T 90: 14-16 s

BOLA Temperature Probes PT 100

» BOLA Immersion Probes PT 100

BOLA Cat. No.	P 1780-20
Materials	Stainless steel tube 1.4571 encapsulated with PTFE, PFA-coated cable
Temperature range	max. +250°C, min. -50°C
Resistance	100 Ohm at 0 °C
Probe-Ø	4 mm
Probe length	approx. 50 mm
Cable length	4.000 mm
Response time	T 50: 4-6 s; T 90: 6-8 s

» BOLA Immersion Probes PT 100 LEMO® either with socket size 1 or plug size 1

BOLA Cat. No.	P 1790-20, P 1792-20
Materials	Stainless steel tube 1.4571 encapsulated with PTFE, PFA-coated cable Chromed brass (LEMO® coupling)
Temperature range	max. +250°C, min. -50°C
Resistance	100 Ohm at 0 °C
Probe-Ø	4 mm
Probe length	approx. 50 mm
Cable length	4.000 mm
Response time	T 50: 4-6 s; T 90: 6-8 s

» BOLA Temperature Probes PT 100 LEMO® Compact with socket size 1

BOLA Cat. No.	P 1730-...
Materials	Stainless steel tube 1.4571 encapsulated with PTFE, Chromed brass (LEMO® coupling)
Temperature range	max. +250°C, min. -50°C
Resistance	100 Ohm at 0 °C
Probe-Ø	8 mm
Probe length	100/300/500 mm
Total length, approx	170/370/570 mm
Probe tip-Ø	5,5 mm
Collar ring-Ø	ca. 12 mm
Response time	T 50: 7-12 s; T 90: 14-16 s

» BOLA Double Temperature Probes PT 100 LEMO® Compact with socket size 1

BOLA Cat. No.	P 1740-...
Materials	Stainless steel tube 1.4571 encapsulated with PTFE, Chromed brass (LEMO® coupling)
Temperature range	max. +250°C, min. -50°C
Resistance	100 Ohm at 0°C
Probe-Ø	8 mm (P 1740-40: 6 mm)
Probe length	300/400/500 mm
Total length, approx	4.000/500/600 mm
Probe tip -Ø	6,0 mm
Collar ring-Ø	27 mm
Response time	T 50: 7-12 s; T 90: 14-16 s

9. Spare Parts / Accessories

» see page 26-27

10. Customized Temperature Probes

Should you be in need of a temperature or immersion probe with other dimensions or connectors, simply send us an inquiry with the requested specification. As manufacturers we are the best contact for custom production.

Generally, the following modifications of our standard temperature probes can be offered:

- » Special dimensions of the thermo element (different probe diameter, useful or total length).
- » Different cable lengths.
- » Different kind of LEMO® connectors – a connector of your choice can be mounted.
- » We can coat your temperature probes with a PTFE heat shrinkable tubing so that they have the chemical resistance of PTFE. Because of the thin coating, the probe or thermometer has slower response behaviour. For coating, your probe has to withstand a short-time temperature of minimum +250 °C.

Please do not hesitate to contact us if you have any questions regarding custom manufacture.

» Registered Trademarks

LEMO® is a registered trademark by Lemo.

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CE Declaration of Conformity

EC-Declaration of Conformity in compliance with Machinery Directive 2006/42/EC
The manufacturer / distributor
Bohlender GmbH · Waltersberg 8 · 97947 Grünsfeld · Germany

herewith declares that the following products:

Temperature Probes Lemo Compact, product group: P 1730- ...

Double Temperature Probes Lemo Compact, product group: P 1740- ...

Temperature Probes, product group: P 1750- ...

Temperature Probes Lemo, product group: P 1760- ... / P 1762- ...

Total Immersion Probes, product group: P 1780- ...

Total Immersion Probes Lemo, product group: P 1790- ... / P 1792- ...

Temperature Probes PT 1000, product group: P 1950- ...

Temperature Probes PT 1000 Lemo, product group: P 1962- ...

Total Immersion Probes PT 1000, product group: P 1980- ...

Total Immersion Probes PT 1000 Lemo, product group: P 1992- ...

apply to the above mentioned directive including any amendments valid at the time this declaration was signed.

The following directives have been applied:

Directive 2014/30 EU

Directive relating to electromagnetic compatibility,
General requirements

The low voltage directive 2006/95/EC has been fulfilled according to 2006/42/EC terms of protections objectives.

The following harmonized standards have been applied:

DIN EN ISO 12100

Safety of machinery

ISO 12100:2010

Risk assessment and risk reduction

Grünsfeld, 03.09.2015

Volker Bohlender

Volker Bohlender

SPECIAL REQUIREMENTS? CUSTOMIZED!

You are looking for something very special? Something that even our huge portfolio of sophisticated lab solutions does not cover?

No problem:

As developer and producer, we offer the possibility to produce individually according to your requirement. This is faster, simpler and often more economic than you can imagine. Just talk to our experts about your ideas – we advise you and support you already during the construction and produce suitable for the material exactly according to your specification. And this already from quantity 1.

For this, we just need a drawing (a rough sketch is sufficient) and some information.

Checklist for your customised product:

- » What is the article name?
- » In which application should the article be used?
- » What dimensions should the article have?
- » Are there any specific material specifications?
- » In which temperature range should the article be used?
- » What chemical stresses is the article exposed to?
- » In which quantities is the article required?
- » What cost per piece should the article not exceed?

You have a special request?
Call us on: **+49 (0) 9346-9286-0**

Or just send us a drawing (a rough sketch is sufficient) and some information to info@bola.de. We will then contact you to discuss further details and establish an offer for you free of charge.

Precise and reliable measurements even in aggressive liquids – all probes are encapsulated with PTFE for maximum chemical resistance.

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BOHLENDER GmbH