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Operating Instructions

Mode d'emploi

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» Thermoelectric Flanging Tools

» Appareils Chauffants

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BOLA Thermoelectric Flanging Tools, Cat.-No: F 701-12 and F 701-14

OPERATING INSTRUCTIONS

1. Supplied with	9
2. Security Advice.....	9
3. Taking into operation.....	9
4. Mode of function.....	11
5. Service / Maintenance.....	11
6. Disturbances / Repairs.....	11
7. Waste disposal.....	13
8. Specification.....	13
9. Spare parts / Accessories	13


Please read the information listed below thoroughly. We urgently ask you to observe the information about security, use and maintenance.

BOLA Thermoelectric Flanging Tools are used to make flanges on tubing made of PTFE, PFA or FEP.

1. Supplied with

- 1 Basic unit (1)
- 1 Flanging tip (assembled) (2)
- 1 Tubing holder (3)
- 1 Operating Instructions

2. Security Advice

- » Do not take a damaged unit into operation.
- » Observe the voltage of the mains: The units with cat.-no. F 701-12 and F 701-14 are suitable only for 230 V / 50 Hz.
- » Unplug the unit for cleaning and maintenance.
- » Do not touch the hot flange tip! Even if the unit has been switched off, there is still a possible risk of burns.
- » Do not use any inflammable cleaning agents!
- » Do not spray any liquids (cleaning agents) directly on the flanging tool.
- » Do not cover the ventilation slots.
- » Do not put any objects into the ventilation slots.
- » Do not work in a moisture area.
- » Risk of injuries due to sharp, pointed flange tip.

3. Taking into operation

- » Check whether the flange tip (2) fits with the inner diameter of the tubing that you will flange.
Flange tips are available for the following tubing I. D.s: 0.5 mm, 0.8 mm, 1.6 mm, 2.4 mm and 3.2 mm. The I. D. is engraved on the top of the flange tip.
- » Exchange only cold flange tips (2).
- » For exchange of the flange tip (2), remove the screw (4) by turning it 2 turns to the left. Use a screw driver blade size 3.5 mm.
- » Draw the flange tip (2) off.
- » Put the new flange tip (2) on the mounting device and push it as far as it will go. Fix it carefully with the screw, otherwise the mounting device can be damaged.

- » Connect the unit with the mains.
Turn the unit on by pressing the switch (7).
The lighting control lamp (7) indicates that the flange tip (2) is heated. The operating temperature is attained after approx. 10 minutes.

Attention! Risk of burns due to hot flange tip (2)!

- » Cut square the tubing to the requested length. Use e. g. the BOLA Tube Cutter (cat.-no. S 1852-28).
Add 2 mm to the total requested length for each flange.

- » Put the Tube End Fitting (5) and the Washer (6) onto the tubing.
The thread of the Tube End Fitting (5) points to the tubing end.

- » Open the included tubing holder (3) by turning the screws. The tube holder is suitable for tubing with outer diameter 1.6 mm (1/16") and 3.2 mm (1/8").


Attention! Risk of burns due to hot flange tip (2)!

- » Place the tubing holder (3) in front of the cooling plate (8). Push the tubing through the tubing holder and as far as it will go. The bores on the cooling plate are for tubing O. D. 1.6 and 3.2 mm. The depth of each bore defines the tubing excess length needed for the flange. Fix the tubing in the tubing holder.

Attention! Risk of burns due to hot flange tip (2)!

- » Push the tubing onto the hot flange tip (2). After approx. 20 seconds the tubing is warm enough and can be formed by applying slight pressure and turning the tube on the flange tip. The exact time has to be defined by making several trials.

- » The flange has to be cooled down fast. Press it for approx. 10 seconds onto the plain surface of the cooling plate (8) (not into the bores). The exact time has also to be defined in several trials.

- » If also the other end of your tubing shall be flanged, do not forget to mount the tube end fitting (5) and the washer (6) first. Then repeat the steps described above.

4. Mode of function

The heating element inside the flanging tool heats the flanging tip until the necessary operating temperature is attained.

A tubing that is pushed onto the hot flange tip gets squashed and the flange can be formed. The flange is cooled down on the cooling plate. The tubing holder facilitates the forming of the flange and helps to press the tube onto the cooling plate.

5. Service / Maintenance

Use a soft and lint-free cleaning cloth. Use a neutral cleaning detergent (ph 7) for strong dirt.

6. Disturbances / Repairs

The following information shall help you to recognise and, if possible, to eliminate possible reasons for disturbances by yourself. If you are not able to handle the disturbance, please give us a call. We will then try to help you this way. If then, a solution to your problem cannot be found, we will ask you to return the cleaned unit in the original packing to us. Then we will check whether a cost-effective repair is possible.

Objection

Reason and help

The tubing end cannot be formed to a flange.

- » Check whether the unit is connected to the mains supply / switched on.
Connect the unit with the mains and switch it on. The control lamp in the switch glows.
- » The assembled flange tip does not fit with your tubing I. D.
Mount the fitting flange tip (the respective size is engraved on the tip).
- » The flange tip has not yet attained the operating temperature.
Heat the flange tip for approx. 10 minutes.
- » The tubing is not warm enough.
The tubing has to be warmed for approx. 20 seconds.
Make several trials to define the exact period of time.
- » Your tubing is not suitable for thermal deforming.
Choose another, suitable tubing.
- » The contact pressure on the flange tip is too low. Press the tubing more strongly to form the flange.


The warmed-up tubing breaks off.

- » The excess length of the tubing end is too long.
Readjust the tubing position in the tubing holder.
- » The warm-up time was too long. Cut the warm-up time short.
- » The tube has not been pressed precisely on the flange tip or on the cooling plate. Press the tubing vertically (right angle) on the flange tip and the cooling plate.
- » The assembled flange tip is too small and thus the tube is not guided during warm-up. Mount a flange tip in the correct size.


The control lamp does not glow although the unit is switched on.

- » The control lamp is faulty.
Ask an electrician to replace the lamp. .
- » The micro-fuse of the power supply is faulty.
The micro-fuse can be exchanged as follows:


- » Pull the plug.

- » Open the fuse holder with a screw driver (blade size 5.5 mm) by turning it counterclockwise.

- » Exchange the faulty fuse against a new micro-fuse type 5/20 0.1 A.
- » Close the fuse holder by turning it clockwise.

7. Waste disposal

BOLA Thermoelectric Flanging Tools as well as their packing are made of recyclable materials. Please take further information on the materials used from the following specifications:

8. Specification

BOLA Thermoelectric Flanging Tools

Cat.-No.	For tubing I. D. mm	Voltage	Tubing holder for tubing O. D. mm	Materials	Power cable length
F 701-12	0.8 mm	230 V / 50 Hz	1.6 mm (1/16") and 3.2 mm (1/8")	Stainless steel, Brass, Aluminium, Hard fibres, Electronic components	2 m
F 701-14	1.6 mm				

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BOHLENDER

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CE Declaration of Conformity

EC-Declaration of Conformity in compliance with Machinery Directive 2006/42/EC

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herewith declares that the following products:

Thermoelectric Flanging Tools, item number : F 701- 12

Thermoelectric Flanging Tools, item number : F 701-14

apply to the above mentioned directive including any amendments valid at the time this declaration was signed.

The low voltage directive 2006/95/EC has been fulfilled according to 2006/42/EC terms of protection objectives.

The following harmonized standards have been applied:


DIN EN ISO 12100

Safety of machinery

ISO 12100:2010

Risk assessment and risk reduction

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Bohlender GmbH • Waltersberg 8 • 97947 Grünsfeld • Germany
BOHLENDER GmbH