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Operating Instructions

Mode d'emploi

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» Flow Filters

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BOLA Flow Filters, Cat.-No.: N 1670-08, N 1670-16 and N 1670-24

OPERATING INSTRUCTIONS

1. SUPPLIED WITH	8
2. SECURITY ADVICE	8
3. FUNCTIONALITY	8
4. TAKING INTO OPERATION	8
5. SERVICE / MAINTENANCE	10
6. DISTURBANCES / REPAIRS	11
7. WASTE DISPOSAL	12
8. SPECIFICATIONS	12
9. ACCESSORIES / SPARE PARTS	12

Please read the information listed below thoroughly. We urgently ask you to observe the instructions about security, use and maintenance.

BOLA Flow Filters can be integrated in an existing line system by means of the included BOLA Laboratory Screw Joints. You can choose the filtering membranes yourself so that they will fit for your own individual requirements. The filtering membranes can easily be assembled or exchanged.

Depending on the version BOLA Flow Filters can be connected to tubing with an outer diameter of 3,2 mm, 6,0 mm, 8,0 mm or 10,0 mm.

BOLA Flow Filters

1. Supplied with

1 BOLA Flow Filter

2 Replacement Caps

4 sets Replacement Inner Parts, consisting of a v-ring, a tapered ring and a sealing ring. Depending on the version the replacement inner parts are for tubing with an O.D. of 3,2 mm, 6,0 mm, 8,0 mm or 10,0 mm.

1 Operating Instructions

2. Security advice


- » Do not take a damaged flow filter into operation.
- » The medium shall not exceed a temperature of $-200\text{ }^{\circ}\text{C}$ up to $+160\text{ }^{\circ}\text{C}$.
- » Please pay attention to the flow direction. A directional arrow is engraved on the housing.
- » The flow rate is reduced if the filtering membrane is saturated.
- » In case the filtering membrane is saturated the pressure of the line system prior to the flow filter can rise respectively the pressure behind the flow filter can decrease.
- » Before integrating the flow filter in your line system make sure that there is no pressure on your line system.
- » Avoid any unwanted reactions with residues if you want to disassemble or open your flow filter.
- » To avoid any health damages we recommend to purge your flow filter with warm water before it is disassembled or before you open it.
- » In accordance to your application please wear suitable protective clothing and safety goggles.

3. Functionality

BOLA Flow Filters can be integrated in an existing line system by means of included Replacement Caps and included Replacement Inner Parts. The optionally available filtering membranes, cat.-no. N 1690-... and N 1564-10, filter any suspended materials out of your medium. The filtrate is only exposed to PTFE.


4. Taking into operation

- » Disassemble the pressure lid (1) by turning it counter-clockwise. Hold the sleeve nut (2).
- » Carefully lift the sealing ring (3) off the filter support (4), e.g. by using a blunt tool like tweezers.
- » **Attention:** The PTFE-coat of the sealing ring must not be damaged.

BOLA Flow Filters


» It's the users liability to choose the right filtering membrane (5) for his application. Put the filtering membrane (5) (available as accessory cat.-no. N 1690-.. and N 1564-10) in the filter support (4). The filtering membrane should lay smooth on the filter support.


» Press the sealing ring (3) into the filter support (4).

» Screw the pressure lid (1) tight into the sleeve nut (2) by turning it clockwise. There is no tool necessary.

» Loosen both replacement caps by one full counter-clockwise turn.


» **Tip:** Sharpen the tubing with a simple sharpener or cut it diagonally

» Slide the tubing through the replacement cap and the inner parts until it bumps against the inside of the flow filter.


» **Tip:** Due to the production tolerances the outer diameter of the tubing might vary. As a result the tubing sometimes can not be slid through the replacement cap and the replacement inner parts. If that's the case you need to slide the components on the tubing one by one. Please pay attention to the order and the position of the replacement cap (6), the v-ring (7), the tapered ring (8) and the sealing ring (9).

BOLA Flow Filters

» Tighten the replacement cap by carefully screwing it counter-clockwise.

» Assemble the other side according to the above mentioned steps.

Flow rate:

Flow capacity under vacuum of 100 kPa (1000 mbar) using a PTFE filtering membrane with a thickness of 0,2 mm:

Membrane-Ø mm	Pore size µm	Product	Flow ml/min
25	1,00	water	96
25	1,00	air	12.800
25	0,45	water	11
25	0,45	air	4.600
47	1,00	water	212
47	1,00	air	48.000
47	0,45	water	50
47	0,45	air	28.000
90	1,00	water	648
90	1,00	air	56.400
90	0,45	water	264
90	0,45	air	36.000

5. Service / Maintenance

Purge the flow filter after usage with a suitable cleaning solvent, e.g. with warm water. For purging the flow filter connect it against the flow direction. If you want to disassemble the flow filter for cleaning proceed as follows:

» **Attention:** Make sure that neither the supply line nor the drainage are under pressure.

» Unscrew both replacement caps by turning them counter-clockwise. Hold the pressure lid (1) and the sleeve nut (2).


BOLA Flow Filters

» The replacement inner parts and the replacement caps may remain on the tubing. This way they can easily be reassembled after cleaning the flow filter.

» Disassemble the pressure lid (1) by turning it counter-clockwise. Hold the sleeve nut (2).


» By using a blunt object, e.g. a pencil, you can push the sealing ring (3), the filtering membrane (5) as well as the filter supporting disc (10) out of the flow filter.

» Reassemble the flow filter after cleansing according to the adjoining exploded drawing.

» You do not need any tools for the assembly. Pay attention to the right mounting direction of the filter supporting disc (10). The filter membrane has to sit on the ring-shaped notches.

6. Disturbances / Repairs

The following information shall help you to recognise and, if possible, to eliminate possible reasons for disturbances by yourself. If you are not able to handle the disturbance, please give us a call. We will then try to help you this way. If then, a solution to your problem cannot be found, we will ask you to return the cleaned flow filter in its original packing to us. Then we will check whether a cost-effective repair is possible.

Objection	Reason and Help
The flow capacity is too small.	<ul style="list-style-type: none"> » The filtering membrane is saturated. Exchange the filtering membrane. » The used filtering membrane has a pore size which is too small for your application. Replace the filtering membrane with a membrane that has a suitable pore size. » The flow direction of the flow filter was overlooked. Thus the filter surface is reduced and therefore less product can flow through the filter. Check the flow direction which is marked by an engraved arrow. » The pressure difference is too low. Increase the pressure prior to the flow filter or reduce the pressure behind the flow filter.

BOLA Flow Filters

The flow filter shows no filter effect.	<ul style="list-style-type: none"> » There is no filtering membrane assembled. Put a suitable filtering membrane into the flow filter. » The used filtering membrane has a pore size which is too big for your application. Use a filtering membrane with a smaller pore size. » The used filtering membrane is damaged. Exchange the filtering membrane with a new one.
Product leaks through the tubing connections.	<ul style="list-style-type: none"> » The tubing is constricted. Disassemble the screw cap. Follow the instructions in point 5 and remove the replacement inner parts. Cut the tubing with a sharp knife for about 5 cm. Attention: Risk of injury! Assemble the tubing connection as described in point 4. » The replacement inner parts have been assembled in a wrong order. Check the order and the position of the v-ring (7), the tapered ring (8) and the sealing ring (9) as described in point 4. » The laboratory screw joint is damaged. Exchange the laboratory screw joint (available as accessory, cat.-no. D 628-..., D 629-... or D 630-...).
Product leaks through the flow filter.	<ul style="list-style-type: none"> » The flow filter is not closed properly. Screw the pressure lid (1) tightly into the sleeve nut (2) by turning it clockwise. There are no tools necessary. » The flow filter cannot be closed completely due to disposals. Purge and disassemble the flow filter as described in point 5. Clean each component individually by using a suitable cleaning solvent, e.g. warm water. Reassemble the flow filter as described in point 4.
The flow filter cannot be opened.	<ul style="list-style-type: none"> » Disposals agglutinate the flow filter. Dive the flow filter with the sleeve nut (2) ahead for about one minute into hot water. Make sure not to dive the pressure lid (1) into the water. Open the flow filter as described in point 5. We recommend to wear heat-proof gloves. Attention: Risk of burns!

7. Waste disposal

The BOLA Flow Filter as well as their packing are made of recyclable materials. Please take further information on the materials used from the following specifications.

8. Specifications

Cat.-No.	Filter-surface	For membrane dia.	Connecting thread	For tubing O.D.	Materials	Temperature resistance
N 1670-08	3,1 cm ²	25 mm	GL 14	3,2 and 6,0 mm	PTFE, Glass-fibre reinforced PPS, FPM (Viton)	- 200 °C up to + 160 °C
N 1670-16	13,8 cm ²	47 mm	GL 18	6,0 and 8,0 mm		
N 1670-24	52,0 cm ²	90 mm	GL 25	8,0 and 10,0 mm		

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BOHLENDER GmbH