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Operating Instructions

Mode d'emploi

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OPERATING INSTRUCTIONS

1. SUPPLIED WITH	6
2. SECURITY INSTRUCTIONS	6
3. TAKING INTO OPERATION	6
4. MODE OF FUNCTION	7
5. SERVICE / MAINTENANCE	7
6. DISTURBANCES / REPAIRS	7
7. WASTE DISPOSAL	8
8. SPECIFICATIONS	8
9. SPARE PARTS / ACCESSORIES	8
10. CUSTOMIZED TEMPERATURE PROBES	8

BOLA Double Temperature Probes
LEMO® Compact ESD

BOLA Temperature Probes PT 100
LEMO® Compact ESD

Please read the information listed below thoroughly. We urgently ask you to observe the information about security and use.

BOLA Double Temperature Probes LEMO® Compact ESD and BOLA Temperature Probes LEMO® Compact ESD are PT 100 probes. They are ideal for measuring temperatures in aggressive liquids and avoid any electrostatic charge. They are available with LEMO® plugs which are mounted directly on the sensor. The plugs are already mounted and can be connected to every suitable measuring instrument by using an extension cable (available as accessory, cat.-no. P 1724-...).

BOLA Temperature Probes Compact ESD

1. Supplied with

1 Double Temperature Probe LEMO® Compact ESD or 1 Temperature Probe PT 100 LEMO® Compact ESD
1 Operating Instructions


2. Security instructions

These temperature probes are designed to make measurements in aggressive liquids and to prevent the build up of an electrostatic charge at the same time. The operating instructions cannot show all security problems which can arise eventually. It's the users liability to make sure that all safety and health regulations are kept and that the corresponding restrictions are stipulated before use.

- » Please do not take a damaged temperature probe into operation.
- » Please take the response times of the temperature probe into consideration during use. Typical response times T 50 and T 90 (50% respectively 90% of the final measuring value) are to be found under point 8 Specifications on page 7.
- » Please do not exceed the operating temperature range of maximum +250 °C and minimum –200 °C.
- » Please note the measuring range of –50 °C up to +250 °C
- » Please verify that the LEMO®-plug is connected correctly to the measuring instrument.
- » Link your earthing holes with your earth circuit connector by means of grounding cable.
- » Please verify the right connection of your grounding cable after each break.

3. Taking into operation

- » Please use an adequate cable (available as accessory, cat-no. P 1724-...) to connect the probe with the measuring instrument.
- » Connect the grounding cable of the Double Temperature Probe LEMO® Compact ESD with an earth circuit connector by using an adequate connecting cable. We recommend a cable cross section of at least 2,5 mm². For attaching the cable please use a screw with a thread diameter of 5 mm and a thread length of 20 mm together with a suitable spring washer and nut.
- » For BOLA Temperature Probes PT 100 LEMO® Compact ESD you can use for example an earth terminal to connect your temperature probe with an earth circuit connector. We recommend a cable cross section of at least 2,5 mm².

BOLA Temperature Probes Compact ESD


4. Mode of function

BOLA Double Temperature Probes LEMO® Compact ESD as well as BOLA Temperature Probes PT 100 LEMO® Compact ESD are sensors based on platinum resistance changes under temperature influence. They are manufactured in conformance to class A. Thus, depending on the real temperature T, the allowed temperature tolerance ΔT is as follows: the deviation of the measured temperature is equal or less than $\pm 0,15 \text{ }^\circ\text{C} + (0,002 \times T)$.

The sensors themselves are located at the end of a stainless steel tube (material code 1.4571). This tube provides certain rigidity, but can be bent to the requested form by hand. The seamless encapsulation with a PTFE/Coal-compound provides an almost universal chemical resistance.

» BOLA Double Temperature Probes LEMO® Compact ESD

This temperature probe has a solid coating made of a conductive PTFE/Coal-compound. By means of a bore below the LEMO® plugs the temperature probe can be connected to an earth circuit connector and thus the build up of an electrostatic charge can be avoided. The surface resistance is 10^4 Ohm . This probe is equipped with two independent, physically separated sensors (1) and (2) (see image on the right). They are manufactured in conformance to DIN EN 61010-2-010:2013.

» BOLA Double Temperature Probes LEMO® Compact ESD

This temperature probe has a solid coating of a conductive PTFE/Coal-compound. By connecting for example an earthing terminal below the LEMO® plug the temperature probe can be connected to an earth circuit connector. Thus the build up of an electrostatic charge can be avoided. The surface resistance is 10^4 Ohm . This probe is equipped with one single sensor.

5. Service / Maintenance

In case of a strong contamination, the temperature probe should be cleaned after use. The PTFE/Coal-compound generally has a smooth, non-wetting surface and can usually be cleaned without any problems. Please do not use abrasive scouring agents because they might damage the surface. You may use all neutral detergents (pH 7). For a stronger contamination we recommend the use of an alkaline detergent up to pH 12. Apart from that, the probes are maintenance-free.

6. Disturbances / Repairs

If the probe should not work properly, please check first if the plug is connected correctly to the measuring instrument. Also a damage of the connecting cable or the thermo element can lead to disturbances. Should there be a disturbance you cannot handle, please do not hesitate to call us. Of course we will try to help you this way. If we do not find a solution to your problem then, we will ask you to return the cleaned probe in the original packing to us. We will check if a cost-effective repair is possible.

BOLA Temperature Probes Compact ESD

7. Waste disposal

The probes as well as their packing are made of recyclable materials. Please take further information about the materials used from the following specifications

8. Specifications

Cat.-No.	Usable length mm	Total length mm ca.	Probe-Ø mm	Number of sensors pieces	Response times	Material	Surface resistance
P 1744-20	300	400	8	2, physically separated, in conformance to DIN EN 61010-2-010:2013	T 50: 20-24 s	PTFE/Coal- Compound	10 ⁴ Ohm
P 1744-23	400	500			T 90: 30 s		
P 1744-30	500	600					
P 1734-10	200	370	1		T 50: 7-12 s		
P 1734-15	300	470			T 90: 14-16 s		
P 1734-20	400	570					
P 1734-25	500	670					

9. Spare Parts / Accessories

» see page 12

10. Customized Temperature Probes

Should you be in need of a BOLA Temperature Probe LEMO® Compact ESD or a BOLA Temperature Probe PT 100 LEMO® Compact ESD with other dimensions or connector couplings, simply send us an inquiry with the requested specification. As manufacturers we are the best contact for custom production.

Generally, the following modifications of our standard temperature probes can be offered:

- » Special dimensions of the thermo element (different probe diameter, useful or total length).
- » Different cable lengths.
- » Different kind of plugs – a plug of your choice can be mounted.
- » We can coat your temperature probes with a PTFE heat shrinkable tubing so that they have the chemical resistance of PTFE. Because of the thin coating, the probe or thermometer has slower response behaviour. For coating, your probe has to withstand a short-time temperature of minimum +250 °C.

Please do not hesitate to contact us if you have any questions regarding custom manufacture.

» Registered Trademarks

LEMO® is a registered trademark by Lemo.

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CE Declaration of Conformity**EC-Declaration of Conformity in compliance with Machinery Directive 2006/42/EC****The manufacturer / distributor**

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herewith declares that the following products:

Temperature Probes Lemo Compact ESD, product group: P 1734- ...

Double Temperature Probes Lemo Compact ESD, product group: P 1744- ...

apply to the above mentioned directive including any amendments valid at the time this declaration was signed.

The following directives have been applied:

Directive 2014/30 EU

**Directive relating to electromagnetic compatibility,
General requirements**

The low voltage directive 2006/95/EC has been fulfilled according to 2006/42/EC terms of protection objectives.

The following harmonized standards have been applied:

DIN EN ISO 12100

Safety of machinery

ISO 12100:2010

Risk assessment and risk reduction

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Bohlender GmbH • Waltersberg 8 • 97947 Grünsfeld • Germany
BOHLENDER GmbH