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Operating Instructions

Mode d'emploi

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» siehe Seite 18

BOLA Pressure-Relief Valve

BOLA Pressure-Relief Valve with manual ventilation, Cat.-No.: E 683-18

Operating Instructions

1. Supplied with	2
2. Security Instructions	2
3. Taking into operation	2
4. Mode of function	3
5. Service / Maintenance	4
6. Disturbances / Repair	5
7. Waste Disposal	6
8. Specifications	6
9. Accessories	6

Please read the information listed below thoroughly. We urgently ask you to observe the information and advice about security, use and maintenance.

BOLA Pressure-Relief Valve

1. Supplied with

- 1 Pressure-relief valve with manual ventilation
- 1 Operating Instructions


2. Security Instructions

- » Do not take a damaged valve into operation.
- » The valve is suitable only for gas or watery products without any solids.
- » Observe the working pressure range from 0 up to 10 bar during operation.
- » The valve is only suitable for product temperatures from $-20\text{ }^{\circ}\text{C}$ up to $+150\text{ }^{\circ}\text{C}$.
- » The flow direction has to be observed. It is marked with an arrow on the housing.
- » If used as an overflow-valve, a controlled deviation of the overflowing medium has to be ensured.
- » Make sure that all components of your installation withstand the pressure that has been adjusted on the valve.
- » Make sure that the requested holding pressure is adjusted.
- » Prior to installation or removing the valve, the inlet and outlet to/from the valve have to be pressureless.


3. Taking into operation

The opening pressure is factory set to 1.5 bar. Any other opening pressure can be adjusted as follows:


- » Loosen the locking nut (1) by turning it counterclockwise upwards the adjusting screw (2).

- » Turn the adjusting screw (2) clockwise with two fingers until you feel a resistance. The pressure-relief valve is now closed.

- » Connect the inlet of the pressure relief valve with a tubing to a compressed air service unit or a pressure gauge. The flow direction is marked with an engraved arrow.
- » Connect the outlet with a tubing (30). Put the tubing end in a vessel with water.
- » Fix the tubes with BOLA Laboratory Screw Joints GL 18.


BOLA Pressure-Relief Valve

» Load the requested opening pressure to the pressure relief valve. It can range between 0.1 and 10 bar.

» Turn the adjusting screw (2) slowly counterclockwise. Watch the water. As soon as bubbles can be seen in the water, the valve is open. The valve's opening pressure is now set to the same pressure as your compressed air service unit.


» Fix the adjusting screw (2) with the locking nut (1). Tighten the locking nut (1) by turning it clockwise with two fingers.

» Check whether the manual ventilation works. Pull the toggle (4). The passage is open. Upon releasing the toggle, the passage is closed immediately.

» **Please note:** the higher the pressure is set, the more power is needed to pull the toggle (4). It can be moved only for a few tenths of a millimetre.

4. Mode of function

The valve is sealed by a membrane with an integrated O-ring. The opening pressure is set by turning the adjusting screw (2) and thus changing the tension applied to the integrated spring. The more the spring is bent respectively the stronger the adjusting screw (2) is closed, the higher is the opening pressure.


By pulling the toggle (4), the membrane is lifted. The passage is open for a short time.

BOLA Pressure-Relief Valve

5. Service / Maintenance

Rinse the pressure-relief valve with clear water after use. Connect the valve to a tap for rinsing. It is possible that the adjusted opening pressure has to be reduced as it has to be lower than the water pressure in the water line. For cleaning of the single components, please observe the following steps for correct disassembly:


» Loosen the locking nut (1). Remove the adjusting screw (2) from the top section (4) by turning it counterclockwise.


» Remove the toggle (3) from the top section (4) by turning it counterclockwise.


» Remove all four screws with a screw driver (blade with 3.5 mm) by turning them counterclockwise.

» Remove the holding plate (5).


» Remove the top section (4).


BOLA Pressure-Relief Valve

- » Remove the pressure plate (6) and the pressure spring (7).
Pull the membrane (8) out of the housing (9).
Tooling is not necessary.

- » Attention: All inner surfaces of the housing (9) and the lower side of the membrane (8) are sealing surfaces. If the surface is damaged, the functionality of the pressure relief valve cannot be guaranteed! Do not remove the O-ring from the membrane.

All single components can be cleaned with commercial, non-abrasive detergents.
Assembly is made in reverse order.

6. Disturbances / Repair

The following information shall help you to recognise and, if possible, to eliminate possible reasons for disturbances by yourself. If you are not able to handle the disturbance, please give us a call. We will then try to help you this way. If then, a solution to your problem cannot be found, we will ask you to return the cleaned valve in the original packing to us. Then we will check whether a cost-effective repair is possible.

A certain minimum pressure shall be maintained in your installation.

Objection	Reason and help
The valve does not open	<ul style="list-style-type: none">» The flow direction of the valve has not been observed. Check the connections and change their order if necessary.» The adjusted opening pressure is too high. It has to be lower than the pressure within your installation. Reduce the opening pressure of the valve.
The valve does not close.	<ul style="list-style-type: none">» There are solids inside the valve. Solids inside the valve can be avoided by adding a filter to the inlet.» The product is too viscous. The pressure relief valve is not suitable for such a product.» The membrane is defective. Exchange it.

BOLA Pressure-Relief Valve

7. Waste Disposal

The BOLA pressure relief valve as well as its packing consist of recyclable material.

Please take further information on the materials used from the following specifications.

8. Specifications

BOLA Pressure Relief Valve with manual Ventilation

Cat.No.	Connecting Thread	Bore dia.	Material	Temperature resistance	Pressure range	Factory setting
E 683-18	GL 18	Ø 6 mm	PTFE PPS stainless steel	from -20 °C to + 150 °C	0,1 bis 10 bar	1,5 bar

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BOHLENDER GmbH