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Operating Instructions

Mode d'emploi

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
BOLA Magnetic Stirrer Head, (G-MRK), Cat.-No.: C 512-08, C 502-08, C 502-16

OPERATING INSTRUCTIONS

1. Supplied with	8
2. Security Instructions	8
3. Mode of Function	8
4. Taking into operation	9
5. Service / Maintenance	10
6. Disturbances / Repairs	11
7. Waste Disposal	11
8. Specifications	12
10. Spare Parts / Accessories	12

Please read the information listed below thoroughly. We urgently ask you to observe the information about security and use.

BOLA Magnetic Stirrer Heads are suitable for stirring in closed vessels with a centre ground joint NS 19, NS 29 and NS 45. The stirrer heads can accept glass and stainless steel stirrer shafts with diameters of 6, 8 and 10 mm. Stirrer shafts with a diameter of 8 and 10 mm must have three counter bores for three stud screws.

The medium is only exposed to metal-free materials which all have a very good chemical resistance.


BOLA Magnetic Stirrer Heads

1. Supplied with

- 1 Operating instruction
- 1 BOLA Magnetic Stirrer Head (G-MRK)


Components of the magnetic stirrer head:

- Pos. 1 Driving bell
- Pos. 2 Conductor
- Pos. 3 Rotor
- Pos. 4 Lower Bearing

2. Security Instructions

- » Do not take a damaged magnetic stirrer head into operation.
- » Only use the magnetic stirrer head up to the maximum speed of 1.500 rpm.
- » Only use the magnetic stirrer head with a stirrer drive which turns to the right (Viewing direction from above: clockwise)
- » Only wear tight clothes to avoid being drawn in by rotating components..
- » If you have long hair, please wear a hairnet.
- » If you use gloves, please make sure that those are tight.
- » Only use the magnetic stirrer head for products which are in a temperature range between -200 and +250°C.
- » Use sleeves with ribs, e.g. from BOLA , Cat.-No.: H 930-...or H 933-.. or grease when you have an application with a reactor lid made of glass.
- » Never use the magnetic stirrer head (B) without globe stirrer coupling (A), Cat. No. C 399-12, or a comparable coupling (see picture on the right).
- » There must be a gap of approx. 2 mm (see picture on the right) between magnetic stirrer head (B) and globe stirrer coupling (A) for compensating heat expansion.
- » Never use the magnetic stirrer heads without the locking ring (5) and the wear rings (6) (see picture on the right).

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4. Taking into operation

Attention:

- » Stirrer shafts with a diameter of 8 mm or 10 mm have to be modified before assembly as shown in the picture on the right side.
- » Stirrer shafts with a diameter of 6 mm must not be modified as they are fixed with a clamping nut.

- » Slide the conductor (2) into the driving bell (1). Make sure that the head of the conductor snaps into the ball bearing of the driving bell properly. The air inside the driving bell may be compressed which can slightly complicate the assembly. Let the air pass off.







- » Slide the lower bearing (4) over your stirrer shaft (5).

Attention: If the stirring diameter of your stirrer shaft (5) is bigger than the ground joint, the shaft (5) has to be pushed from below through the reactor lid before the lower bearing (4) and the rotor can be assembled.

- » For stirrer shafts with diameter 8 or 10 mm:
Remove the three stud screws (7) in the rotor (3) by turning them counter-clockwise. Use a screwdriver with a blade tip of 4 mm.


- » Place the counter bores of the stirrer shaft (5) in parallel to the threaded holes of the rotor (3). Screw the stud screws (7) into the counter bores of the shaft (5) by turning them clockwise. Make sure that the stud screws (7) are completely screwed into the rotor (3).

- » For stirrer shafts with diameter 6 mm:
Open the clamping nut (8) by turning it counter-clockwise for at least two times.


BOLA Magnetic Stirrer Heads

- » Push the stirrer shaft up against the rotor (3).
- » Close the clamping nut (8) tightly by turning it clockwise with two fingers.


- » Push the stirrer shaft with the lower bearing (4) and the assembled rotor (3) from below into the conductor (2). Be aware that the rotor (3) is sucked into the conductor on the last few centimetres due to the magnetic pull.

- » Tighten the lower bearing (4) in the conductor (2) by turning it counter-clockwise.

Attention: left-handed thread!

- » etzen Sie Ihre Magnet-Rührkupplung vorsichtig in Ihren Schliff.
- Tip:** For female ground joints made of glass, we recommend to use sleeves with ribs, e. g. BOLA Cat.-No. H 930-... or an equivalent. If your application allows it, you may use grease.

Specifics of the BOLA Magnetic Stirrer Head

- » Using the magnetic stirrer head at high speed can cause heat development on the conductor (2) and on the driving bell (1). The components will get hand-hot.
- » In non-operating state, there is a tiny space between the lower bearing (4) and the rotor (3). This space may close during use.
- » If the magnetic stirrer head is waved a rattling sound can be heard. Those specifics do not affect the operation.

5. Service / Maintenance


Please do not use abrasive detergents since these might roughen the surface and therefore damage it. You can use all common neutral detergents (pH 7).

For a stronger contamination we recommend to use an alkaline detergent up to pH 12.

The BOLA Magnetic Stirrer Heads can be used without any additional lubricants

BOLA Magnetic Stirrer Heads

» During permanent use at high speed particles of abrasion can stuck on the level of the wear rings (6) inside the conductor. These can be removed after disassembling the magnetic stirrer head and has no affect on the operation. For cleaning use only non-abrasive detergents.

» The wear rings (6) on the rotor (3) and on the lower bearing (4) as well as the locking ring (5) should be visually inspected every time the stirrer head is being disassembled. A deformed locking ring (5) has to be exchanged, as it does not secure the lower bearing anymore. Consequently the lower bearing (4) might get lose while the stirrer head is in use. The wear rings (6) have to be changed if they are thinner than 1,0 mm.

6. Disturbances / Repairs

The following information shall help you to recognise and, if possible, to eliminate possible reasons for disturbances by yourself. If you are not able to handle the disturbance, please give us a call. We will then try to help you this way. If then, a solution to your problem cannot be found, we will ask you to return the cleaned magnetic stirrer head in the original packing to us. Then we will check whether a cost-effective repair is possible.

Objection	Reason and Help
The driving bell (1) is moving, the rotor (3) stands still.	» The locking ring (5) on the lower bearing (4) was not assembled, therefore the rotor (3) is clamped. Assemble the locking ring (5).
The transmitted torque is significantly lower than the given one.	» The driving bell (1) has not the right position to the rotor (3). Hence the magnets cannot transmit the full force. Check the position of the driving bell(1) on the conductor (2). Make sure that the head of the conductor (2) snaps into the ball bearing of the driving bell (1) properly.

7. Waste Disposal

The BOLA Magnetic Stirrer Head as well as its packing consist of recyclable material. Please take further information on the materials used from the following specifications.

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CE Declaration of Conformity

EC-Declaration of Conformity in compliance with Machinery Directive 2006/42/EC

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herewith declares that the following products:

Magnetic Stirrer Heads (G-MRK), product group: C 502- ...

Magnetic Stirrer Heads (G-MRK), product group: C 504- ...

Magnetic Stirrer Heads (G-MRK), product group: C 512- ...

apply to the above mentioned directive including any amendments valid at the time this declaration was signed.

The following directives have been applied:

Directive 2014/30 EU

**Directive relating to electromagnetic compatibility,
General requirements**

The low voltage directive 2006/95/EC has been fulfilled according to 2006/42/EC terms of protection objectives.

The following harmonized standards have been applied:


DIN EN ISO 12100

Safety of machinery

ISO 12100:2010

Risk assessment and risk reduction

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Bohlender GmbH • Waltersberg 8 • 97947 Grünsfeld • Germany
BOHLENDER GmbH