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Operating Instructions

Mode d'emploi

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- » siehe Seite 13

BOLA Hydrolyzing and Digestion Vessels

BOLA Hydrolyzing and Digestion Vessels, Cat. No.: A 240-02, A 240-04, A 240-06, A 240-08 and A 240-10.

BOLA Hydrolyzing and Digestion Vessels with Liners, Cat. No. A 250-04, A 250-06 and A 250-08.

Operating Instructions

1. Supplied with	6
2. Security instructions	6
3. Taking into operation	6
4. Mode of function	7
5. Waste disposal	7
6. Service / Maintenance	7
7. Disturbances / Repairs	8
8. Specifications	8
9. Spare parts /Accessories	8

Please read the information listed below thoroughly. We urgently ask you to observe the information about security, use and maintenance.

BOLA-Hydrolyzing and Digestion Vessels are made to be used in microwave ovens. They are made of the fluoroplastic PTFE TFM 1600. Its homogeneous, non-porous and extremely smooth surface avoids so-called memory-effects. Due to the compact structure as well as the almost universal chemical resistance they are especially well suited for digestions with aggressive mediums with a pressure up to 25 bars (depending on the volumes).

BOLA Hydrolyzing and Digestion Vessels

1. Supplied with



- 1 Hydrolyzing and Digestion Vessel made of PTFE TFM 1600
- 1 Liner (only included with Cat. No. A 250-04, A 250-06 or A 250-08)
- 10 Sealing membranes made of PFA
- 10 Rupture membranes made of PTFE
- 1 Operating instructions

2. Security instructions

- » Please do not take a damaged vessel into operation.
- » Never use the Digestion Vessel with Liner without liner.
- » Never use the digestion vessel without collecting vessel.
- » Please observe the values for pressure and temperature (see specifications).
- » Please consider that the material of the vessel can be warmed by the content.
- » The digestion vessel should only be opened after a sufficient cooling period of at least 20 minutes.
- » Place the digestion vessel under a fume hood for opening.
- » Wear protective clothing and goggles when opening the digestion vessel.
- » Please consider that the content of the vessel can still be hot even if the outer surface is cold.
- » Please renew sealing and rupture membrane for every digestion.
- » Never use metalliferous liquids for digestion in the microwave.
- » Never heat the digestion vessel on a heating plate or with a flame.
- » Please observe the respective security instruction when working with acids and solvents.

3. Taking into operation

- » Unscrew the lid (1) from the vessel (2) by turning counterclockwise.
- » If you have a digestion vessel with liner, please check if the liner (10) is assembled. The liner can be pushed into the vessel (11) with little pressure.
- » Push the rupture membrane (8) made of PTFE (milky) into the lid (1).


BOLA Hydrolyzing and Digestion Vessels


- » Fill the vessel (2) with your sample.
- » Place a sealing membrane (9) made of PFA (clear) on the vessel (2). Close the lid (1) tightly by turning clockwise.
- » Release the nut (3) by turning counterclockwise (approx. 2 turns). The screw joint is suitable for a tubing O.D. of 6.35 mm (1/4").

- » Push up the tubing (4) against the nut (3).
- » Close the nut of the screw joint (3) by turning clockwise with two fingers.
- » Put the end of the tubing into a collecting vessel (for example BOLA Jars, Cat. No. A 131-12).

- » **Tip:** The O.D. of a tubing can vary due to production tolerance. If the tubing cannot be pushed into the screw joint, we recommend to assemble the components separately. Please observe the arrangement of sealing ring (6) and tapered ring (7).

4. Mode of function

The vessel is closed tightly with the massive lid. The inserted sealing membrane seals up between vessel and lid. The inserted rupture membrane is acting as "pressure relief valve" and protects the digestion vessel from overpressure. The screw joint with assembled tubing leads liquids into a collecting vessel. In a microwave oven, only the content of the digestion vessel is heated.

5. Waste disposal


The packing as well as the vessels are made of recyclable materials.

6. Service / Maintenance

Please use commercial, non-aggressive detergents for cleaning. An intensive cleaning of the TFM parts can be made by steaming them out first with hydrochloric or nitric acid and then with water. Please consider following instruction for removing and cleaning the liner:

BOLA Hydrolyzing and Digestion Vessels

- » Push the liner (10) with a blunt item (20) out of the housing (11), then pull it out completely.

7. Disturbances / Repairs

The following information shall help you to find by yourself the reason and a possible solution for disturbances. In case of disturbances you cannot handle, please do not hesitate to contact us by phone. Of course we will try to help you and give advice to find a solution. If the problem might not be solved this way, the cleaned and originally packed vessel can be shipped to us. We will check then if a cost-effective repair is possible.

Problem	Reason and help
Vessel leaks	<ul style="list-style-type: none"> » Check the sealing membrane and change it if necessary. » Upper and lower part of the vessel are not screwed up enough, please screw up stronger. » The sealing surface between housing and sealing membrane is damaged.
Sealing and rupture membranes burst.	<ul style="list-style-type: none"> » The pressure was too high. Check the microwave oven concerning the time regulated and the power. Check the quantity of the content and reduce it if necessary.

8. Specifications

BOLA Hydrolyzing and Digestion Vessels

Cat. No.	Material	Capacity	Max. pressure	Max. temperature	Internal dimensions	External dimensions
A 240-02	PTFE TFM	5 ml	25 bar	160 °C	Ø 15 x 32 mm	Ø 40 x 93 mm
A 240-04	1600	10 ml			Ø 16 x 52 mm	
A 240-06	PTFE	20 ml	20 bar	150 °C	Ø 22 x 60 mm	Ø 95 x 159 mm
A 240-08	PFA	50 ml			Ø 33 x 62 mm	
A 240-10		100 ml	15 bar	140 °C	Ø 35 x 110 mm	Ø 95 x 191 mm

BOLA Hydrolyzing and Digestion Vessels with Liners

Cat. No.	Material	Capacity	Max. pressure	Max. temperature	Internal dimensions	External dimensions
A 250-04	PTFE TFM	10 ml	25 bar	160 °C	Ø 24 x 63 mm	Ø 60 x 125 mm
A 250-06	1600	20 ml			Ø 30 x 63 mm	
A 250-08	PTFE	50 ml	20 bar	150 °C	Ø 43 x 77 mm	Ø 95 x 159 mm

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BOHLENDER GmbH