




Sophisticated Safety Solutions for UHPLC.





b.safe Caps UHPLC



bsafe.de/en

b.safe Caps UHPLC

Through their freely movable insert, they can be easily screwed without disarranging the connected capillaries.

Article	Scope of delivery	Cat.No.
Cap UHPLC	1 x b.safe Cap 1 x Air Valve 3 x Fitting 3,2 mm 2 x Blind Fittings	M 1045-03



Chemically highly resistant

With the cap made of robust PPS and the chemically highly resistant PTFE insert, b.safe Caps GL45 can be sterilized and autoclaved and easily cleaned in the dishwasher.




Best grip

The ergonomic design and eroded surface of the b.safe Cap UHPLC also provides a firm grip for gloved hands of all sizes.

b.safe Air Valves UHPLC



bsafe.de/en



b.safe Air Valves UHPLC

The upstream membrane with a pore size of 0.1 μm made of porous PTFE effectively retains coarse and fine dust, as well as microorganisms and other impurities carried by it from the outside.

Article	Scope of delivery	Cat.No.
Air Valve UHPLC	2 x Air Valve	M 551-01

- PTFE filter membrane (pore size 0.1 μm)
- No way for microorganisms to get through
- No algae and biofilm growth
- No blockage of capillaries and filter frits
- No downtimes, no maintenance costs

b.safe

Where I can get more informations?

bsafe.de

+49 (0) 9346 9286-0

info@bohlender.de

f



in



ig



Follow us