

## Strong Novelties

Besides our well-proven classics you will find again many innovations in the BOLA catalogue. On the following pages we show you a choice of our novelties. Detailed information can be found on the stated pages.



### **BOLA Globe Stirrer Couplings**

**Cat. No.: C 394-..**

For balancing misalignment  
of axes between agitator and  
stirrer shaft



### **BOLA Vacuum Adaptors with Ground Joint**

**Cat. No.: N 1655-..**

To put on standard  
filter funnels

# Continuation - Novelties



## **BOLA Vacuum Adaptor GL**

Cat. No.: N 1656-..

Adaptor made of PTFE to put on standard filter funnels

Page 323



## **BOLA Ground Joint Connector**

Cat. No.: B 306-..

Connection of two components with ground joint sockets

Page 212



## **BOLA Multiple Distributors with Ground Joint**

Cat. No.: D 623-..

Easy extension of connection options

Page 216



## **BOLA T-Threaded Couplings**

Cat. No.: H 908-..

Connection of three GL threads with the same thread size

Page 127



## **BOLA Reactor Lids with Center Ground Joint Neck**

Cat. No.: B 120-..

Lid for all standard reaction vessels with flat flange

Page 267



## **BOLA Flat Flange Lids**

Cat. No.: B 210-..

Lid for all standard reaction vessels with flat flange

Page 280



## **BOLA PTFE Tubing Reels**

Cat. No.: S 2010-.. / S 2012-..

PTFE tubing ready-made in different lengths

Page 193



## **BOLA FEP Tubing Reels**

Cat. No.: S 2020-.. / S 2022-..

FEP tubing ready-made in different lengths

Page 194



## **BOLA PFA Tubing Reels**

Cat. No.: S 2030-.. / S 2032-..

PFA tubing ready-made in different lengths

Page 194



## **BOLA PTFE-EX Tubing Reels**

Cat. No.: S 2040-..

PTFE tubing ready-made in different lengths

Page 157