

b.safe SYSTEM

Professional equipment
for fill level control

b.safe Mini-Display



Better b.safe!

With the new **b.safe Mini-Display**.

✓ optical
fill level control

✓ easy to operate
via sensor surface

✓ low installation
depth

✓ acoustic
fill level control

✓ suitable for installation
in control panels

✓ flexible connection
cables

”

With the b.safe Mini-Display you are always informed in time with visual and acoustic signals. This means that containers no longer overflow unnoticed. In addition, they are compatible with all b.safe fill-level sensors.

“



Advantages of the **b.safe Mini-Display**

Optical signal

To visually indicate the fill level, the LED lights up in different colours



Acoustic signal

As soon as the maximum fill level is reached, an acoustic signal sounds



Low installation depth

Installation is possible from an installation depth of 32 mm thanks to the extended connection cable.



Minimum wall thickness for installation

Can be mounted on wall thicknesses from 1 to 35 mm. Suitable for panel mounting.

Easy to operate via touch sensor surface

With snooze function.



Precise fill level sensors

The b.safe Mini-Display is compatible with all b.safe fill level sensors.



Simple installation

No laborious plugging in thanks to the flexible and extended cables. The extended connectors allow easy connection to the device even under difficult conditions.

b.safe SYSTEM

Our product selection



b.safe Mini-Display Wall
1-channel
Cat. No.
R 658-01



b.safe Mini-Display
Table 1-channel
Cat. No.
R 659-01



b.safe Mini-Display Wall
2-channel
Cat. No.
R 658-02



b.safe Mini-Display
Table 2-channel
Cat. No.
R 659-02



b.safe Mini-Display
DLC Wall 1-channel
Cat. No.
R 658-03



b.safe Mini-Display
DLC Table 1-channel
Cat. No.
R 659-03



b.safe Mini-Display
Splitter Wall, 2-channel
Cat. No.
R 658-04



b.safe Mini-Display
Splitter Table, 2-channel
Cat. No.
R 659-04

Where else can I get information?

bsafe.de

+49 (0) 9346 9286-0

info@bohlender.de



Follow us