

# b.safe SYSTEM

Professional equipment  
for fill level control

## b.safe Mini-Display



# Better b.safe!

With the new **b.safe Mini-Display**.

✓ optical  
fill level control

✓ easy to operate  
via sensor surface

✓ low installation  
depth

✓ acoustic  
fill level control

✓ suitable for installation  
in control panels

✓ flexible connection  
cables

”

With the b.safe Mini-Display you are always informed in time with visual and acoustic signals. This means that containers no longer overflow unnoticed. In addition, they are compatible with all b.safe fill-level sensors.

“



# Advantages of the **b.safe Mini-Display**


## Optical signal

To visually indicate the fill level, the LED lights up in different colours



## Acoustic signal

As soon as the maximum fill level is reached, an acoustic signal sounds



## Low installation depth

Installation is possible from an installation depth of 32 mm thanks to the extended connection cable.



## Minimum wall thickness for installation

Can be mounted on wall thicknesses from 1 to 35 mm. Suitable for panel mounting.

## Easy to operate via touch sensor surface

With snooze function.



## Precise fill level sensors

The b.safe Mini-Display is compatible with all b.safe fill level sensors.



## Simple installation

No laborious plugging in thanks to the flexible and extended cables. The extended connectors allow easy connection to the device even under difficult conditions.

# b.safe SYSTEM

## Our product selection



b.safe Mini-Display Wall  
1-channel  
**Cat. No.**  
**R 658-01**



b.safe Mini-Display  
Table 1-channel  
**Cat. No.**  
**R 659-01**



b.safe Mini-Display Wall  
2-channel  
**Cat. No.**  
**R 658-02**



b.safe Mini-Display  
Table 2-channel  
**Cat. No.**  
**R 659-02**



b.safe Mini-Display  
DLC Wall 1-channel  
**Cat. No.**  
**R 658-03**



b.safe Mini-Display  
DLC Table 1-channel  
**Cat. No.**  
**R 659-03**



b.safe Mini-Display  
Splitter Wall, 2-channel  
**Cat. No.**  
**R 658-04**



b.safe Mini-Display  
Splitter Table, 2-channel  
**Cat. No.**  
**R 659-04**

Where else can I get information?

[bsafe.de](http://bsafe.de)

+49 (0) 9346 9286-0

[info@bohlender.de](mailto:info@bohlender.de)



Follow us